
KACLIP®

GBRUCHSANWEISUNG • INSTRUCTIONS FOR USE

DE KANÜLENTRAGEBAND	HU KANÜLTARTÓ SZALAG
EN TUBE HOLDER	PL TASMA STABILIZACJA RURKIĘ
FR BANDE DE FIXATION DE CANULE	RO LEINTA DĂRĂ FIXACĂIJA KANIOIUI
IT FASCIA DI FISSAGGIO CANNULA	CS PŮLNĚN KANÝLY
ES CINTA DE FIJACIÓN DE LA CÁNULA	SK PÁS NA NOSIENIE KANÝLY
PT FAIXA DE FIXAÇÃO DA CÁNULA	SL TRAK ZA PROČISTITEV KANILE
NL CANULEDRAGBAND	HR FIKSACIJA VRPCE ZA KANULU
SV KANÜLYBAND	SR POVEZNA TRAKA ZA KANULU
DA KANYLEBÆREBÆND	BG ЗАКРЕПВАЩА ЛЕНТА ЗА КАНОЛИАТА
NO BÆRESTROPP FOR KANYLE	RO BANȚA DE SUSȚINERE A CANULEI
FI KANNYLINNOIKKOHINA	JA カニューライト
EL ΤΑΜΙΑ ΣΤΕΡΕΩΣΗΣ ΤΡΑΧΕΙΩΣΗΣ/ΗΝΗ	AR شريط حمل الكانيل
TR KANÜLÜ YÖNÜ BAĞI	

BILDER / PICTURES

KACLIP®

Place the tube holder around your neck with the hooks/hook-and-loop fasteners facing outward **⊗** ab. Attach the two hooks/hook-and-loop fasteners to the lateral eyelets on both sides of the neck flange **⊗** ab. Adjust the slide on the tubeholder in accordance with your individual neck size **⊗**.

VI. PERIOD OF USE
The tube holder is a single-patient product and also suitable for repeated use.
The service life depends on the individual frequency of use as well as the condition and on individual handling of the product.
Place the tube holder in use must be placed temporarily on a clean surface.
Before being reattached, the product must be checked for visible damage and for functionality.
The product may not be reused if it is damaged, heavily soiled or heavily soaked, as this could jeopardize safe use.
In this case, the product must be replaced immediately.
The tube holder must be disposed of and replaced with a new tube holder no later than 7 days after first legal notice.

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EN I. FOREWORD
These instructions for use are valid for tube holders. The instructions for use serve to inform the physician, nursing staff and patient/user in order to ensure technically correct handling. Please read the instructions for use carefully before using the product for the first time!

II. INTENDED USE
The tube holder allows the tracheostomy tube to be fixed securely in place and can be adjusted to fit the individual neck size of the patient.

III. WARNINGS
The tube holder is a single-patient product and is thus intended exclusively for use - including repeated use - in one and the same patient. The tubeholder should not be fitted too tightly around the neck. The slide should be at a distance of approx. 1 cm between the tube holder and the neck. To prevent the tracheostomy tube from coming off by accident, the tube holder must be checked regularly to make sure it is still securely fixed in place.

IV. INSTRUCTIONS
Place the tube holder around your neck with the hooks/hook-and-loop fasteners facing outward **⊗** ab. Attach the two hooks/hook-and-loop fasteners to the lateral eyelets on both sides of the neck flange **⊗** ab. Adjust the slide on the tubeholder in accordance with your individual neck size **⊗**.

V. PERIOD OF USE
The tube holder is a single-patient product and also suitable for repeated use.
The service life depends on the individual frequency of use as well as the condition and on individual handling of the product.
Place the tube holder in use must be placed temporarily on a clean surface.
Before being reattached, the product must be checked for visible damage and for functionality.
The product may not be reused if it is damaged, heavily soiled or heavily soaked, as this could jeopardize safe use.
In this case, the product must be replaced immediately.
The tube holder must be disposed of and replaced with a new tube holder no later than 7 days after first legal notice.

VI. LEGAL NOTICES
The manufacturer Andreas Fahí Medizintechnik-Vertrieb GmbH assumes no liability for damage caused by unauthorised product modifications, repairs not carried out by the manufacturer or by improper handling, care (cleaning/disinfection) and/or storage of the products contrary to the provisions of these instructions for use. Should a manufacturer's liability in connection with this product of Andreas Fahí Medizintechnik-Vertrieb GmbH, this is to be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient are domiciled. The manufacturer reserves the right to make changes to the product without notice. KACLIP® is a trademark of Andreas Fahí Medizintechnik-Vertrieb GmbH, Cologne, which is registered in Germany and the European Member States.

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VI. LEGAL NOTICES
The manufacturer Andreas Fahí Medizintechnik-Vertrieb GmbH assumes no liability for damage caused by unauthorised product modifications, repairs not carried out by the manufacturer or by improper handling, care (cleaning/disinfection) and/or storage of the products contrary to the provisions of these instructions for use. Should a manufacturer's liability in connection with this product of Andreas Fahí Medizintechnik-