








TRACHEOSILC® / TRACHEOTEC® SILC

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PICTROGRAM LEGEND

EN

| | | | |
|--|--|---|-------------------------|
|  | Outer cannula with no inner cannula |  | Single patient use only |
|  | Outer cannula including 15 mm-swivel connector (VARIO) |  | Do not reuse |
|  | With suction device (SUCTION) |  | Keep dry |
|  | CUFF |  | Do not resterilize |
|  | MRT suitable |  | Medical devices |
|  | Obturator (insertion aid) | | |
|  | Tube holder included | | |
|  | BE For ventilation | | |
|  | Use by | | |
|  | Date of manufacture | | |
|  | Manufacturer | | |
|  | Content (in pieces) | | |
|  | Consult instructions for use | | |
|  | Do not use if package is damaged | | |
|  | LOT Batch code | | |
|  | STERILE EO Sterilized using ethylene oxide | | |
|  | REF Catalogue number | | |

TRACHEOSILC® & TRACHEOTEC® SILC TRACHEOSTOMY TUBES

I. FOREWORD

These instructions for use are valid for Fah!® tracheostomy tubes. The instructions for use serve to inform the physician, nursing staff and patient/user in order to ensure technically correct handling of Fah!® tracheostomy tubes.

Please read the instructions for use carefully before using the product for the first time!

Keep them in an easily accessible place for future reference.

Please keep the package for as long as you use the tracheostomy tube. It contains important information on the product!

II. INTENDED USE

Fah!® tracheostomy tubes are intended to stabilise the tracheostoma following laryngectomy or tracheostomy.

The tracheostomy tube is designed to keep the tracheostoma open.

Tracheostomy tubes with cuff are always indicated after tracheostomy has been performed for any cause if the cannula needs to be sealed against the wall of the trachea.

When used for the first time the products must only be selected, used and inserted by clinicians or medical professionals trained in their use.

III. WARNINGS

Patients must be trained by medical professionals on how to use the Fah!® tracheostomy tubes safely.

Fah!® tracheostomy tubes must never be blocked, e.g. by mucus or encrustation. Blockage can cause asphyxiation!

Mucus in the trachea can be aspirated (removed by suction) through the tracheostomy tube with a tracheal suction catheter.

Damaged tracheostomy tubes must not be used and must be discarded immediately. Use of a damaged tube may result in airway compromise.

When the tracheostomy tubes are inserted or removed, this can lead to irritations, coughing or slight bleeding. If bleeding persists, consult your doctor immediately!

Do not use tracheostomy tubes during laser treatment (laser therapy) or with electrosurgical devices! Contact with a laser beam may damage the tube.

Cleaning, disinfection, re-sterilisation and re-use can impair the safety and function of the product and are therefore not permissible!

CAUTION!

Tracheostomy tubes that contain metal parts must under no circumstances be used during a radiation therapy (radiotherapy), as this can lead to serious skin damage! If you need to wear a tracheostomy tube during a radiation therapy, then only use plastic tracheostomy tubes without any metal parts. In the case of plastic tracheostomy tubes with silver speaking valve, the valve can for instance be removed completely from the tube, including the safety chain, by removing the inner cannula from the outer cannula before starting a radiation therapy.

IV. COMPLICATIONS

The following complications could emerge when using this device:

Contamination of the stoma may make it necessary to remove the tube. Contaminations can furthermore cause infections that necessitate the use of antibiotics.

Unintentional Inhalation of a tube that was not sized correctly may need to be removed by a doctor/physician. If mucus blocks the tube it should be removed and cleaned.

V. CONTRAINDICATIONS

Do not use if the patient is allergic to the material.

CAUTION!

Never under any circumstances use cannula variants without cuff for patients on mechanical ventilation!

VI. PRECAUTIONS

The correct tracheostomy tube size should be selected by the treating physician or by trained medical professionals.

When using Fah!® tracheostomy tubes, the products may become infected for instance with yeast (Candida), bacteria or the like, which may lead to the material becoming cloudy and its service life being shortened. Should this occur, the tube must be replaced immediately.

In order to rule out accidental detachment of auxiliary equipment or damage to the tube, insert only auxiliary equipment with a 15 mm connection into Fah!® UNI adapters.

We highly recommend always keeping at least two replacement tubes on hand to ensure continuous supply.

CAUTION!

When increased tensile forces occur during mechanical ventilation, e.g. due to stiff swivel connectors connected to the tracheostomy tube or due to uncontrolled movements of the patient, the inner cannula can be accidentally twisted out of the outer cannula. The patient must therefore be monitored, or the tracheostomy tube and/or the tube system must be replaced if necessary.

Do not perform any type of adjustment, repair or alteration on a tracheostomy tube or compatible accessories. If a product is damaged, it must be discarded immediately.

VII. PRODUCT DESCRIPTION

The TRACHEOTEK[®] SILC and TRACHEOSILC[®] tracheostomy tubes are products which are made of medical-grade silicone.

The tracheostomy tubes are made of thermosensitive medical synthetics that develop their optimal product properties at body temperature.

We supply the Fahl[®] tracheostomy tubes in various sizes and lengths.

Please see the size table in the appendix.

TRACHEOTEK[®] SILC tracheostomy tubes have a high-pressure cuff. Please heed in this regard Section 3.1.2 on leakproofness testing of the cannula and high-pressure cuff.

TRACHEOSILC[®] / TRACHEOTEK[®] SILC Tracheostomy tubes are single-patient products and are intended for single use.

The package contains 1 tube, supplied sterile by ethylene oxide sterilization (EO).

The tube tip is rounded to prevent irritation of the mucus membranes in the trachea.

It may be advisable to use tubes of different lengths alternately so the tip of the tube does not always touch the same spot in the trachea and cause irritation. This prevents pressure points or the formation of granulation tissue. Always discuss the exact procedure with your doctor/physician.

MRT (Magnetic Resonance Tomography) INFORMATION

These tracheostomy tubes can even be worn by the patient during magnetic resonance tomography (MRT).

This only applies to tracheostomy tubes without metal parts / without cuff.

CAUTION!

Since the tracheostomy tubes with cuff have a small metal spring in the non-return valve of the pilot balloon with inflation tube, the tracheostomy tube with cuff must not be used while an MRI (magnetic resonance imaging) examination is being performed.

MRT is a diagnostic imaging technique in which inner organs, tissues and joints are depicted with the help of magnetic fields and radio waves. Metallic objects can be drawn into the magnetic field and can induce changes by their acceleration. Although the metal spring is extremely small and lightweight, interactions that could give rise to health impairments or to malfunction or damage of the technical equipment used or of the tracheostomy tube itself nevertheless cannot be ruled out in this context. If wearing a tracheostomy tube is indicated for keeping the tracheostoma open, we recommend use of a metal-free tracheostomy tube instead of the tracheostomy tube with cuff while an MRT is being performed.

CAUTION!

Do not use a tracheostomy tube with cuff while a nuclear magnetic imaging / MRI examination is being performed!

1. Neck flange

A particular feature of the Fahl[®] tracheostomy tube is the specially shaped neck flange, which is designed to fit the anatomy of the neck.

The size specifications are printed on the neck flange.

The neck flange of the tracheostomy tubes has two side eyelets for inserting a tube holder.

All Fahl[®] tracheostomy tubes with fastening eyelets include a disposable tube holder. The tube holder secures the tracheostomy tube to the neck.

Please follow the instructions provided with the tube holder carefully when attaching and removing it from your tracheostomy tube.

The insertion aid (obturator) included in delivery facilitates insertion of the tracheostomy tube.

Ensure that the Fahl[®] tracheostomy tube is placed in the tracheostoma without tension and the position is not changed when fastening the tube holder.

2. Connectors and adapters

Connectors and adapters are used to connect compatible tube accessories.

The options available to a patient depend on the clinical condition, such as status post laryngectomy or tracheostomy.

Connectors/adapters are normally permanently attached to the inner cannula. The connector in this case is the universal (15 mm) connector with the aid of which so-called artificial noses or heat and moisture exchangers (HMEs, filters for exchanging heat and moisture) can be attached.

A 15 mm standard connector enables secure connection to required accessories.

3. Cannula tube

The cannula tube borders directly on the neck flange and directs the airflow into the trachea.

The tube tip is rounded to prevent irritation of the mucus membranes in the trachea.

3.1 Cuff

In the product versions with cuff, the very thin-walled, high-volume cuff adapts well to the trachea and ensures reliable sealing if inflated correctly. The cuff can be inflated like a balloon. The small pilot balloon on the inflation tube indicates whether the tracheostomy tube is in sealed (inflated) or non-sealed condition.

The cuff itself is inflated via a tube with one-way valve and pilot balloon.

CAUTION!

Remove the safety protector of the valve on the pilot balloon if such is present (see picture 17).

3.1.1 Leakproofness test of the cannula and the cuff (if such exists)

The cannula and the cuff must be tested for leakages directly before and after every insertion and thereafter at regular intervals.

For this purpose, fill the cuff with 18 to 22 mmHg (1 mmHg corresponds to 1.35951 cmH₂O) and observe whether a spontaneous pressure drop occurs.

There should be no substantial drop in pressure in the cuff during the observation period.

This leakproofness test must also be performed prior to every renewed insertion (e.g. after cleaning the tracheostomy tube) (see picture 7c).

3.1.2 Leakproofness test of the cannula and the high-pressure cuff (if such exists)

The high-pressure cuff is inflated exclusively using a syringe. The optimal inflation of the high-pressure cuff has been reached when no leakage can be heard anymore or the ventilator machine no longer indicates a leakage.

CAUTION!

The inflation volume of the high-pressure cuff must be specified by the physician.

Never under any circumstances use a cuff pressure gauge to inflate the high-pressure cuff.

Observe the maximum inflation volume of the cuff (see table of sizes).

The following are possible signs indicating leakages in the cuff (balloon):

- Externally visible damage to the balloon (holes, fissures etc.)
- Audible and/or tactile perception of air escaping from the balloon
- Water in inlet tubes of the tracheostomy tube (after cleaning!)
- Water in the cuff (after cleaning!)
- Water in the pilot balloon (after cleaning!)
- No cough stimulus when pressure is applied to the pilot balloon

CAUTION!

When testing the balloon or when inserting or removing the tracheostomy tube, never under any circumstances use sharp or pointed objects, such as for instance forceps or clamps, since these can damage or even destroy the balloon. Should one of the above-named signs of leakage be detectable, the tracheostomy tube must under no circumstances be used any more since it is no longer functioning properly!

3.2 Obturator

Please check whether the obturator can easily be removed from the cannula before inserting the tracheostomy tube! After checking the obturator for ease of movement, reinsert the obturator into the cannula for inserting the tracheostomy tube.

The insertion aid serves to stabilise the tracheostomy tube during insertion into the tracheostoma.

Insertion aids are single-patient products and are intended for single use only. They must not be cleaned or disinfected.

3.3 Suction opening (only in Suction tracheostomy tube versions)

Secretions that have accumulated above the inflated cuff can be removed via the suction opening in the outer tube of the Suction tracheostomy tubes.

The suction tube leading to the outside can be connected to a syringe or a suction device. The user must decide in consultation with the treating physician and after appropriate risk analysis, which of these suctioning variants to use. The individual disease condition of the patient should be considered in this context.

Suctioning using a suctioning device may in any case only be performed if the suctioning device is equipped with vacuum adjuster. The suction pressure must never exceed – 0.2 bar.

When using the suction variants of the tracheostomy tubes, the secretions that have accumulated above the cuff must be suctioned off via the suction opening on the tracheostomy tube directly before deflating the cuff in order to prevent aspiration of secretions. In this way, additional suctioning with a suction catheter while deflating the cuff may be rendered unnecessary, which facilitates handling significantly.

CAUTION!

When using any of the tracheostomy tube versions with a suctioning device („SUCTION“), take special care during the suctioning procedure to generate a vacuum for as short a time as possible; drying out of the subglottic space as a result of this must be avoided.

CAUTION!

Contraindicated in patients with increased bleeding tendency (e.g. during anticoagulant therapy). In such cases, the Suction tracheostomy tube with suction opening may not be used, since suctioning entails an increased risk.

VIII. DIRECTIONS FOR TUBE INSERTION AND REMOVAL**For the Doctor/Physician**

The appropriate tracheostomy tube must be selected by a doctor/physician or trained medical professionals.

Select a tube that fits the patient's anatomy to optimise comfort and ventilation (breathing in and out).

For the Patient**CAUTION!**

Always insert the tracheostomy tube with the cuff in fully deflated condition (see picture 7a)!

CAUTION!

Carefully examine the sterile packaging to ensure that it has not been tampered with or damaged. Do not use the product if the packaging has been damaged.

Check the use-by or expiry date. Do not use after this date.

It is advisable to use sterile disposable gloves.

Carefully examine the tube before first use to make sure that it is not damaged and that there are no loose parts.

Should you notice any anomaly or anything unusual, DO NOT use the tube. Return the tube to the manufacturer for inspection.

The neck flange respectively the funnel-shaped housing (retainer ring) must not be inserted into the tracheostoma. Make sure that it is always outside the tracheostoma (see picture 2).

If secretion collects in the lumen of the Fah!® tracheostomy tube and cannot be removed by coughing or aspiration, the tube should be removed.

1. Insertion of the tube

Step-by-step instructions to insert Fah!® tracheostomy tubes.

Before application, users should clean their hands (see picture 3).

Remove tube from the package (see picture 4).

If an obturator is to be used, this must first be fully inserted into the cannula tube so that the collar on the gripping piece of the obturator comes into contact with the outer edge of the 15 mm connector and the tip of the olive projects beyond the tip of the cannula (proximal end of cannula). The obturator must be held in this position during the entire procedure.

When using tracheostomy tubes with cuff, pay special attention to the following points:

Before inserting the tracheostomy tube, check the cuff (balloon) as well - the cuff must be free of damages of any kind and must be leakproof so as to ensure tight sealing as required. We therefore recommend to perform a leakproofness test prior to every insertion (see Section VII, 3.1.1 / 3.1.2).

The balloon must be completely empty prior to insertion of the tracheostomy tube (see picture 7b)! When using an aid for dilating the tracheostoma, take care to ensure that the tracheostomy tube, and especially the cuff, is not damaged by friction.

Then a tracheal compress is pushed onto the tracheostomy tube.

To increase the lubricity of the tracheostomy tube and thus facilitate insertion into the trachea, it is recommended to wipe the outer tube with an OPTIFLUID® stoma oil wipe (REF 31550), which allows even distribution of the stoma oil on the tube (see Figs. 4a and 4b) or FAHL® OPTIFLUID® Lubricant Gel 20g tube (REF 36100) or FAHL® OPTIFLUID® Lubricant Gel 3g sachet (Ref 36105).

If you are inserting the tube yourself, use a mirror to make insertion of the Fah!® tracheostomy tube easier.

When inserting the Fah!® tracheostomy tube, hold it by the neck flange with one hand (see picture 5). Pull the tracheostoma slightly apart with your free hand to allow the tip of the tube to fit into the tracheostoma more easily.

Special aids for dilating the tracheostoma are also available (tracheal dilator, REF 35500). These allow the tracheostoma to be dilated gently and evenly, for instance also in emergency situations with collapsing tracheostoma (see picture 6).

Make sure that the tube is not damaged by friction when using an instrument for assistance.

Now carefully insert the tube into the tracheostoma during the inspiration phase (while breathing in) while tilting your head slightly back (see picture 7).

Advance the tube into the trachea.

Straighten your head once the tube has been inserted further into the trachea.

If an obturator is used, this must then immediately be removed from the tracheostomy tube.

The tracheostomy tubes should always be attached with a special tube holder. This stabilises the tube and thus ensures that the tracheostomy tube is securely seated in the tracheostoma (see picture 1).

1.1 Inflating the cuff (if such exists)

The cuff is inflated by applying a defined pressure to the cuff via the Luer connection (standardised conical connection) of the inflation tube by means of a cuff pressure gauge. If not instructed otherwise by the doctor/physician, we recommend a cuff pressure of at least 18 mmHg (20 cmH₂O) to 22mmHg (25 cmH₂O).

The high-pressure cuff is inflated by applying a defined pressure to the cuff via the Luer connection (standardised conical connection) of the inflation tube by means of a syringe.

Inflate the cuff at most to this target pressure and check to make sure that sufficient air is supplied via the tracheostomy tube.

Always make sure that the cuff is undamaged and is functioning faultlessly.

If the desired sealing is not achieved even after repeated attempts with the specified limit volume, a tracheostomy tube with larger diameter may be indicated

The correct cuff pressure must be checked regularly, i.e. at least every 2 hours.

CAUTION!

All instruments used for inflating the cuff must be clean and free of foreign particles! Detach the instruments from the Luer connection of the inflation tube as soon as the cuff has been inflated.

CAUTION!

If the maximum pressure is exceeded for longer periods of time, the blood circulation in the mucus membrane can be impaired (risk of ischaemic necrosis, pressure ulcers, tracheomalacia, tracheal stenosis, pneumothorax). In patients undergoing artificial respiration, the cuff pressure should not be allowed to drop below the cuff pressure value specified by the doctor/physician in order to prevent unnoticed aspiration. Hissing noises in the region of the balloon, especially during expiration, indicate that the trachea is insufficiently sealed by the balloon. If the trachea cannot be sealed with the pressure values specified by the doctor/physician, the entire air should be withdrawn again from the balloon and the sealing process should be repeated. If this does not lead to success, we recommend to use the next larger tracheostomy tube with balloon. Due to the permeability of the balloon wall for gases, it is normal for the pressure in the balloon to decline slightly over time, but it can on the other hand also rise unintentionally during gas anaesthesia. Regular pressure monitoring is therefore urgently recommended.

The cuff must never under any circumstances be inflated with excessive amounts of air, since this can lead to damage of the tracheal wall, tears in the cuff with subsequent deflation, or deformation of the cuff, in which case airway obstruction cannot be ruled out.

CAUTION!

During anaesthesia, the cuff pressure can rise/fall due to nitrous oxide (laughing gas).

2. Removing the tube

CAUTION!

Accessories such as a tracheostoma valve or HME (Heat Moisture Exchanger) must be removed first before proceeding to remove the Fahl[®] tracheostomy tube.

CAUTION!

If the tracheostoma is unstable, or in emergency situations (puncture/dilation tracheostomy), the tracheostoma can collapse after withdrawal of the tracheostomy tube, thereby impairing air supply. A fresh tracheostomy tube must be kept ready for use in such cases and must be quickly inserted if necessary. A tracheal dilator (REF 35500) can be used for temporarily securing the air supply.

The cuff must be emptied before removing the tracheostomy tube. The head should be tilted back slightly for removal of the tube.

CAUTION!

Never use a cuff pressure gauge to empty the cuff. Always use a syringe for this.

Before the air is removed from the balloon by means of a syringe and the tracheostomy tube is withdrawn, the region of the trachea above the balloon must first be cleaned by suctioning off secretions and mucus. If the patient is responsive and reflexes are intact, it is recommended that the patient be suctioned while at the same time unblocking the tracheostomy tube. Suctioning is performed by inserting a suction catheter through the cannula tube into the trachea. In this way, suctioning can be performed without any problems and gently for the patient and cough stimulus and the risk of aspiration are minimised.

Then deflate the cuff while suctioning at same time. If secretions are present, these are now removed and can no longer be aspirated. Please note that prior to reinsertion, the tracheostomy tube must always be cleaned, disinfected if necessary and lubricated with stoma oil as specified below.

Proceed very carefully to avoid injury to the mucus membranes.

Step-by-step instructions to remove the Fahl[®] tracheostomy tubes:

The tracheostomy tubes should be removed with the head slightly tilted back. Grip the tube at the side by the neck flange or the housing (see picture 7). Carefully remove the tracheostomy tubes.

IX. STORAGE/CARE

Still sterile packed replacement tubes should be stored in a dry environment away from dust, sunlight and heat.

Only use stoma oil (OPTIFLUID® Stoma Oil, 25 ml bottle REF 31525/ Stoma oil wipe REF 31550) or lubricant gel (FAHL® OPTIFLUID® Lubricant Gel 20 g tube REF 36100/ 3 g sachet REF 36105) as lubricant.

We highly recommend keeping at least two replacement tubes on hand to ensure continuous supply.

X. SERVICE LIFE

These tracheostomy tubes are sterile products for single-patient use.

The maximum period of use should not exceed 29 days.

The service life of a tube is influenced by many factors. The composition of the secretions and other aspects, for example, are very important.

Damaged tubes must be replaced immediately.

CAUTION!

Modifications to the tracheostomy tube of any kind, in particular perforations as well as repairs to the tube, may only be carried out by the manufacturer himself or by companies expressly authorised for this in writing by the manufacturer! Inexpertly performed work on tracheostomy tubes can lead to serious injuries.

XI. LEGAL NOTICES

The manufacturer Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for functional deficiencies, injuries, infections, and/or other complications or adverse events caused by unauthorised product alterations or improper use, care, and/or handling.

In particular, Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for damages caused by modifications to the tracheostomy tube, especially due to shortening of the same or perforations, if these modifications or repairs were not carried out by the manufacturer himself. This applies both to damages to the tubes themselves caused thereby and to any consequential damages caused thereby.

If the tracheostomy tube is used beyond the period of use specified under Section X, and/or if the tracheostomy tube is used, maintained or stored in non-compliance with the instructions and specifications laid down in these instructions for use, Andreas Fahl Medizintechnik-Vertrieb GmbH will, as far as legally permissible, be free of any liability, including liability for defects.

Should a serious adverse event occur in connection with this product of Andreas Fahl Medizintechnik-Vertrieb GmbH, this is to be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient are domiciled.

Sale and delivery of all Andreas Fahl Medizintechnik-Vertrieb GmbH products is carried out exclusively in accordance with our General Terms and Conditions of Business which can be obtained directly from Andreas Fahl Medizintechnik-Vertrieb GmbH.

Product specifications described herein are subject to change without notice.

TRACHEOSILC®/TRACHEOTEC® is a trademark and brand of Andreas Fahl Medizintechnik-Vertrieb GmbH, Cologne registered in Germany and the EU member states.

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All measurements listed in the tables below are for standard cannulas.
Customised cannulas are available upon request.

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Sizing information:

All lengths in mm | GR = Size | AD1 = Outer diameter at the neck flange | AD2 = Outer diameter at cannula tip | ID1 = Inner diameter at the cannula tip of the outer cannula | ID2 = Inner diameter at the cannula tip of the inner cannula | L1 AB = Length over the outer curve | L2 MB = Length over the middle curve | Ø-CD = CUFF diameter | Ø-BW = Bending angle

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