Instruction for use

FİLTRELER FILTERS

IFU ID: AW.04.ALT.AL.0868-REV.08 (03/2020)

Contents

AR	تعليمات التشغيل (العربية).....	2/12
CZ	Český (Návod K Použití).....	2/12
DA	Dansk (Brugervejledning).....	3/12
DE	Deutsch (Bedienungsanleitung).....	3/12
EL	Ελληνικά (Οδηγος Χρησης).....	4/12
EN	English (Instructions For Use).....	4/12
ES	Español (Instrucciones De Uso).....	5/12
FA	دستورالعملهای عملیاتی (فارسی).....	5/12
FR	Français (Mode D'emploi).....	6/12
HU	Magyar (Használati Utasítás).....	6/12
IT	Italiano (Istruzioni Per L'uso).....	7/12
NL	Nederlands (Gebruiksaanwijzing).....	7/12
PL	Polski (Instrukcje Użytkowania).....	8/12
RO	Romana (Instrucțiuni De Utilizare).....	8/12
RU	Русский (Инструкция По Применению).....	9/12
TR	Türkçe (Kullanım Talimatları).....	9/12
ZH	中文 (使用说明).....	10/12

Manufactured by:

Meditera Tibbi Malz. San. ve Tic. A.Ş.

İbni Melek OSB Mh. Tosbi Yol 4 Sok. No:29

35900 Tire İzmir

TURKEY

Tel: (+90 232) 513 51 10 (Pbx) www.meditera.com.tr

BRUGSANVISNING FOR FILTER

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ENGLISH

USER MANUAL FOR FILTER

The following instructions are provided for maximum performance and safety in use of breathing filters.

Description:

Single patient use HME and antibacterial/viral filter for breathing circuits in anesthesia, reanimation and intensive care for bacterial and viral filtration as well as humidification of patient airways.

Instruction for use:

Connect the device securely to breathing system, catheter mount, anesthesia/ventilator machine, closed suction system ensuring that the connections are tight and safe. Also check if present perfect seal of the cap for capnometry port.

Warnings:

Please perform control tests for any blockage and leakage on the product and all respiratory line before starting the treatment.

1. The filter has been manufactured for single patient use and maximum recommended hours of use is 24 hours.
2. This filter is only used between the patient and breathing circuit.
3. The product is manufactured under clean room conditions. Do not use if the product and product package is damaged.
4. This product should be used only and exclusively under supervision of a physician.
5. Replace the filter in the event of increased resistance due to secretion from the patient.
6. The additional dead space of this device must be taken into consideration when used. The product is latex-free.

Indication	
	Do not reuse
 eIFU indicator	Consult Ifu.
	Do not use if package is damaged.
	Keep away from sunlight.
	Sterilized using ethylene oxide.
	Do not restrelize
	Keep away from rain
	CE mark. The product complies with the EU Council directive concerning Medical Devices 93/42/EEC.

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CE işareti. Ürün, Medikal Cihazlar ile ilgili AB Konsey direktifi 93/42/EEC'ye uygundur.