
NEBUFIRST® PRO

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
STANDARD ACCESSORIES SUPPLY

- 1 - Nosepiece
- 2 - Mouthpiece
- 3 - Adult mask
- 4 - Pediatric mask
- 5 - Medicine Cup - topside
- 6 - Medicine Cup - underside
- 7 - Medicine Cup - air entrance
- 8 - Pisper
- 9 - Connecting Tube
- 10 - Air filter for replacement

en

PISTON NEBULIZER FOR AEROSOLTHERAPY

SAFETY WARNINGS

THE DEVICE IS A MEDICAL DEVICE FOR HOME USE ONLY.

IT MUST BE PUT INTO OPERATION FOLLOWING TO THE LETTER WHAT IS INDICATED IN THIS MANUAL.

IT IS IMPORTANT THAT THE USER READ AND UNDERSTAND THE INFORMATION FOR THE USE AND MAINTENANCE OF THE UNIT.

CONTACT THE MANUFACTURER WITH ANY QUESTIONS

THE INTENDED USE OF THIS DEVICE IS EXCLUSIVELY FOR AEROSOL THERAPY

- In the presence of children and non-self-sufficient people, the appliance must be used under the supervision of an adult who has read this manual and has full mental faculties.
- In the presence of children and non-self-sufficient people, the appliance must be used under the supervision of an adult who has read this manual and has full mental faculties.
- The ampoule and accessories are supplied NOT sterile. Before use, read the "CLEANING AND MAINTENANCE" paragraph
- The ampoule and accessories must be reused after proper cleaning. Personal use is recommended
- The device is a Class II a Medical Device that transforms drugs or solutions from liquid form to aerosol form, allowing them to be administered by air. Nebulizes drugs both in suspension and in solution
- The therapeutic effect can be reduced with particularly oily and / or viscous substances
- It is not a nebulizer for essential oils
- Do not spray liquids and / or alcohol-based preparations
- The device generates compressed air.
- Unwanted side effects can arise from intolerance to the type of pharmaceutical speciality used.
- Our company policy provides for continuous improvement; therefore the packaged product may be slightly different from the one illustrated
- Carefully read the warnings contained in this manual and keep it carefully for future reference
- After removing the appliance from the packaging, make sure that it is intact without visible damage.
- Otherwise, contact only technically competent personnel or the point of purchase
- Use the device only as a nebulizer for aerosol therapy.
- Any use other than that intended is improper and therefore dangerous for the user.

- In case of improper use, the manufacturer is not to be considered responsible.
- Device not suitable for anesthesia and pulmonary ventilation.
- The correct operation of the device could be compromised by electromagnetic interference due to a malfunction of nearby devices such as television, radio, etc.
- If this happens, try moving away from the offending device until the interference disappears, or try connecting the aerosol to an alternative electrical outlet
- Never use extension cables or adapters.
- It is recommended to always unwind the power cable to avoid dangerous overheating.
- Keep the cord away from hot surfaces
- Always disconnect the plug from the socket after use
- Do not immerse the device in water; should this happen, immediately disconnect the plug from the socket
- Do not remove or touch the appliance immersed in water before having disconnected the plug from the socket.
- Do not reuse the appliance after removing it from the water and immediately contact the technical assistance center
- Do not use the appliance with wet hands
- Do not obstruct the cooling grids in any way
- Make sure that the appliance during use is positioned on flat and stable surfaces to avoid accidental falls and overturns
- Do not use the device while the patient is sleeping
- For any repairs, contact the manufacturer or qualified personnel only.
- Do not open or tamper with the device.
- If the appliance does not work properly, turn it off and consult the instruction manual
- Do not leave the packaging components (plastic bags, cardboard boxes, etc.) within reach of children as they are potential sources of danger
- Keep nebulizer unit and accessories out of the reach of children.
- Children and non-self-sufficient people must use the appliance only under the strict supervision of an adult who has read this manual
- Do not leave the device with the plug inserted in the socket when not in use
- Personal use of accessories is recommended
- Do not leave drug residues in the neblizer unit and / or accessories once the therapy is finished
- Cleaning and / or maintenance must take place only after turning off the appliance and unplugging the power cord from the power outlet
- Do not pull the power cord or the appliance itself to remove the plug from the socket.
- The device is a medical device for home use only.

- It must be put into operation following to the letter what is indicated in this manual.
- It is important for the patient to read and understand the information for the use and maintenance of the unit.
- Contact the manufacturer with any questions
- Do not obstruct the air outlet
- Check the instructions before use:
 - **0 / OFF** Switch off
 - **I / ON** Switch on
 - ~ Alternating current
- Do not use the appliance while bathing or showering
- Store the device in a dry, clean place away from heat sources and atmospheric agents
- Do not use the appliance near sources of potential interference (example: heat sources such as fireplaces, electric or gas stoves; sources of steam such as kettles, sources of electromagnetic radiation such as microwave ovens)
- Keep away from the reach of pets that could contaminate surfaces in contact with the patient, obstruct the ventilation grids, and prevent the correct course of therapy.
- The product contains small parts that can be inhaled or swallowed
- The power cord can be a source of strangulation risk: supervise children and non-self-sufficient patients during therapy
- The device is to be used only with the medicines and the relative dosage indicated by your doctor and the information sheet of the medicine
- The appliance has been designed for domestic use.
- The maximum continuous operation is 20 minutes.
- The minimum pause after operation is 40 minutes
- Variations in the number of revolutions of the motor, due to voltage drops, etc., do not modify the average performance of the device
- Use only the accessories supplied with the product, other accessories can make the therapy ineffective
- Using the device in violation of one or more of the warnings described above could cause the device to malfunction.
- Incorrect intake of the therapy prescribed by the attending physician and damage to the patient
- Full compliance with medical warnings and prescriptions is recommended to avoid serious problems

For the correct disposal of the device, comply with the regulations in force

0051 The device complies the MDD 93/42/CE Directive.

INSTRUCTIONS FOR THE PREPARATION AND START-UP OF THE DEVICE

Not being sterile, the device needs cleaning and disinfection before and after each use.

(Follow the instructions in the paragraph "INSTRUCTIONS FOR CLEANING THE DEVICE").

To insert the solution to be nebulized, unscrew the upper part as indicated in figure "A" and pour the dose indicated by the doctor into the lower part (B).

Make sure that the pisper is inserted in the lower part as indicated in figure "C".

Then screw the upper part back on and connect the tube to the lower part (D), taking care not to tilt the latter to avoid liquid leaks.

Connect the other end of the tube to the aerosol therapy device.

Place the accessory you want to use from those supplied with the (E) nebulizer on the upper part.

Turn on the device connected to the ampoule following the instruction manual supplied with the device to start the therapy.

⚠ ATTENTION!

KEEP THE NEBULIZER UNIT IN A VERTICAL POSITION DURING THE TREATMENT; DO NOT TILT MORE THAN 45°

THE PRODUCT CONTAINS SMALL PARTS THAT MAY BE SWALLOWED OR INHALED: KEEP THE PRODUCT OUT OF THE REACH OF CHILDREN

WARNINGS FOR THE CORRECT CONDUCT OF THERAPY

- Sit comfortably and upright
- If you are using the mask (adult or pediatric), make it adhere well to the face.
- The mask is equipped with special holes for the outflow of air during the exhalation phase.
- Breathe slowly with a short pause at the end of each breath
- If you are using the mouthpiece, hold it lightly between your teeth with your lips closed.
- Breathing must take place through the mouth
- If you are using the nose fork, place it gently on the nostrils, not inserting it inside them. Breathing must take place through the nose
- The therapeutic application is finished when the ampoule emits the nebulized particles in a discontinuous way
- Once the application is complete, carry out a thorough washing and / or disinfection as described in the "CLEANING AND DISINFECTION" chapter



Fig. E

en

CLEANING AND DISINFECTION

After each therapy cycle, unplug the device from the power outlet and clean the device and accessories used:

CLEANING OF THE NEBULIZER UNIT, NOSEPIECE AND MOUTHPIECE

Clean the components of the nebulizer unit, the mouthpiece and the nosepiece using hot water (max 60 °C) for about 5 minutes, possibly adding a minimum quantity of detergent, taking care to remove any residual medicine and dirt deposits (Observe the instructions of the detergent manufacturer).

Rinse thoroughly, remove excess water and leave to dry in a clean place. It is possible to disinfect the nebulizer unit with chemical disinfectants readily available in pharmacies suitable for use on plastic materials (Amuchina®, Milton®, etc.).

Follow the instructions, dosage indications and use restrictions provided by the disinfectant manufacturer.

After disinfection is complete, allow the device to return to room temperature if hot water is used before using it again.

CLEANING THE MASK AND THE CONNECTION TUBE

Thoroughly wash only with chemical disinfectant.

MICROBIAL CONTAMINATION

In the presence of pathologies with risks of infection and microbial contamination, personal use of the accessories and nebulizer kit is recommended.

REPLACING THE FILTER

To replace the filter, lift the filter holder cover by levering it with a flat screwdriver (E).

Remove the worn filter, clean the cover with a damp non-abrasive cloth and insert a new replacement filter among those supplied with the product. Put the cover with the new filter (F) back in its seat by pressing lightly.



USEFUL LIFE OF THE DEVICE

This device has a useful life of approximately 2,500 therapy cycles, considering an average therapy duration of 8 minutes.

The nebulizer unit supplied with the device must be replaced every 100-120 therapy cycles.

The filter should be replaced every 30 hours of use or when it takes on a greyish color.

For information on replacing the filter, refer to the paragraph "REPLACING THE FILTER".

TECHNICAL SPECIFICATIONS

- Rotary piston compressor, lubrication free
- Class IIa medical device according to MDD Directive 93/42 / EEC annex IX
- Linear ampoule particle size: < 5 µm MMAD
- Intermittent use: 20 min ON – 40 min OFF
- Nominal voltage: 230 V ~
- Frequency: 50Hz
- Consumption: 120 VA
- ENVIRONMENTAL CONDITIONS OF USE
- Ambient temperature: min 5 ° C - max 40 ° C
- Air humidity: min 15% RH - max 93% RH
- Atmospheric pressure: min 700hPa - max 1060hPa (maximum operating altitude ≤ 2000m)

en


CONDITIONS OF TRANSPORT AND STORAGE

- Temperature min -25 ° C / max + 70 ° C
- Air humidity: min 0% RH - max 93% RH
- Atmospheric pressure: min 500 hPa - max 1060 hPa

SYMBOLS ADOPTED

In compliance with CEI EN 60601-1: 2006 + A1: 2013, CEI UNI EN ISO 15223-1: 2018, ISO 7001: 2007 and ISO 7000: 2004 standards, the symbols used are the following:

~	AC VOLTAGE		CE MEDICAL MARKINGS REF. DIR. 93/42 EEC
	DOUBLE INSULATION DEVICE		DO NOT USE THE APPLIANCE WHILE BATHING OR SHOWERING
	TYPE BF APPLIED PART	I / ON	SWITCH ON
	OBLIGATION TO READ THE INSTRUCTION MANUAL	O / OFF	SWITCH OFF
	GENERAL WARNING SYMBOL		INDICATES THE MANUFACTURER OF THE MEDICAL DEVICE
IP21	PROTECTION FROM THE PENETRATION OF LIQUID AND SOLID BODIES (PROTECTED AGAINST SOLID BODIES LARGER THAN 12MM; PROTECTED AGAINST ACCESS WITH A FINGER; PROTECTED AGAINST THE VERTICAL FALL OF WATER DROPS)		THIS PRODUCT (INCLUDING ITS ACCESSORIES) AT THE END OF ITS USEFUL LIFE, MUST BE TAKEN TO A SEPARATE COLLECTION CENTER FOR ELECTRICAL AND ELECTRONIC EQUIPMENT BY THE USER, OR BY YOUR TRUSTED DEALER

	INDICATES THE NEED FROM PART OF THE USER OF SEE INSTRUCTIONS FOR USE.		INDICATES THE LOT NUMBER ON THE RATING LABEL
	INDICATES THE REFERENCE CODE OF THE DEVICE ON THE RATING LABEL		INDICATES THE SERIAL NUMBER OF THE PRODUCT ON THE RATING LABEL
	INDICATES THE CORRECT VERTICAL POSITION OF THE TRANSPORT PACKAGING.		INDICATES THE RANGE OF PRESSURE VALUES ATMOSPHERIC TO WHICH THE MEDICAL DEVICE CAN BE EXHIBITED IN SAFETY.
	INDICATES A DEVICE DOCTOR WHO CAN BREAK OR BE DAMAGED IF IT IS NOT HANDLED CAREFULLY.		INDICATES THE LIMITS OF TEMPERATURE AT WHICH THE MEDICAL DEVICE CAN BE EXHIBITED IN SAFETY.
	INDICATES THAT THE EXTERNAL PACKAGING OF THE PRODUCT IS RECYCLABLE		INDICATES A MEDICAL DEVICE WHO NEEDS HUMIDITY PROTECTION.

PROBLEMS AND SOLUTIONS

1. The appliance does not turn on:
 - Make sure that the plug is firmly inserted in the socket
 - Make sure that the switch is in the on position
2. The device is ON but does not nebulize:
 - Make sure you have inserted the pisper inside the cruet
 - Make sure that the air connection hose is not squashed or kinked
 - Make sure the air intake filter is not clogged or dirty.
If it is, replace it with a new one
 - Make sure that the drug is inside the nebulizer unit.

If the appliance does not resume proper operation despite the checks carried out, contact the manufacturer

IMPORTANT INFORMATION CONCERNING ELECTROMAGNETIC COMPATIBILITY

Due to the increasing number of electronic devices (PCs, tablets, smartphones, etc.), the medical devices in use may be subject to electromagnetic interference produced by other equipment.

Such interference could cause the device to malfunction and this could create an unsafe situation.

Furthermore, medical devices must not interfere with other equipment.

For compliance with electromagnetic compatibility (EMC) regulations and in order to prevent potential risk situations in the use of the product, the standards contained in EN60601-1-2: 2015 + A1: 2021 have been implemented.

These standards define the levels of immunity to electromagnetic interference, as well as the maximum levels of electromagnetic emissions for medical devices.

The medical devices produced by EP S.p.A. comply with the standards EN60601-1-2: 2015 + A1: 2021 with regard to both immunity and emissions.

However, the following precautions must be observed:

- Do not use the device near cell phones, smartphones or other devices that generate strong electric or electromagnetic fields.
This could cause the unit to malfunction and create a situation potentially unsafe.
- It is recommended to keep such equipment at a minimum distance of 7 m.
- Check the correct operation of the device if the distance is less

Further documentation relating to the standards EN60601-1-2: 2015 + A1: 2021 is available from EP S.p.A. at the address indicated in this instruction manual.

Electromedical devices require particular care, with regard to EMC requirements, it is therefore required that they be used in accordance with the manufacturer's specifications.

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Andreas Fahl Medizintechnik-Vertrieb GmbH