INHEALTH

TECHNOLOGIES®

InHealth 22.5 Fr Indwelling Voice Prosthesis

INSTRUCTIONS FOR USE

37918-01B Effective January 2020 / 2020 年 1 月生效 / 2020年1月発効 / Në fuqi nga janari 2020 / 2020 اعتبارًا من شهر يناير / В сила от януари 2020 г. / Stupa na snagu u siječnju 2020. / V platnosti od ledna 2020 / Gældende fra januar 2020 / Met ingang van januari 2020 / Kehtib alates jaanuarist 2020 / Voimassa tammikuusta 2020 alkaen / Date d'entrée en vigueur : janvier 2020 / Gültig ab Januar 2020 / Σε ισχύ από τον Ιανουάριο του 2020 / 2020 בתוקף החל מינואר / Hatályos 2020. januártól / Gildir frá janúar 2020 / Valido da gennaio 2020 / 2020년 1월 발효 / Spēkā no 2020. gada janvāra / Įsigalioja 2020 m. sausio mėn. / Gjelder fra januar 2020 / Wchodzi w życie w styczniu 2020 r. / Vigente em janeiro de 2020 / Cu intrare în vigoare din ianuarie 2020 / Версия от января 2020 г. / Važi od januara 2020. / Účinné od januára 2020 / Efectivo a partir de enero de 2020 / Gäller från januari 2020 / Geçerlilik Tarihi: Ocak 2020

InHealth Technologies® is a registered trademark of **Freudenberg Medical, LLC**.

Patent(s): www.inhealthpatents.com







Copyright © 2020 Freudenberg Medical, LLC. All rights reserved.

INSTRUCTIONS FOR USE

使用说明 / 使用說明書 / Udhëzime për përdorim / إرشادات الاستخدام / Инструкции за употреба / Upute za korištenje / Návod k použití / Brugsanvisning / Gebruiksaanwijzing / Kasutusjuhend / Käyttöohjeet / Instructions d'utilisation / Gebrauchsanweisung / Οδηγίες χρήσης / הוראות שימוש / Használati utasítás / Notkunarleiðbeiningar / Istruzioni per l'uso / 사용 지침 / Lietošanas norādījumi / Naudojimo instrukcijos / Bruksanvisning / Instrukcja użytkowania / Instruções de uso / Instrucțiuni de utilizare / Инструкции по применению / Uputstvo za upotrebu / Návod na použitie / Instrucciones de uso / Bruksanvisningar / Kullanım Talimatlar

DIAGRAMS

插图 / 図 / Diagramet / الرسم التوضيحي / Схеми / Dijagrami / Schémata / Diagrammer / Afbeeldingen / Diagrammid / Kaaviot / Diagrammes / Abbildungen / Διαγράμματα / תרשימים / Ábrák / Skýringarmyndir / Diagrammi / 그림 / Diagrammas / Diagramos / Diagrammer / Schematy / Diagramas / Diagrame / Схеми / Crteži / Grafické znázornenia / Diagramas / Diagram / Diyagramlar







TABLE OF CONTENTS

目录 / 目次 / Tabela e përmbajtjes / جدول المحتويات / Съдържание / Sadržaj / Obsah /
 Indholdsfortegnelse / Inhoud / Sisukord / Sisällysluettelo / Table des matières /
 Inhaltsverzeichnis / Πίνακας περιεχομένων / תוכן עניינים / Tartalomjegyzék / Efnisyfirlit /
 Sommario / 목차 / Saturs / Turinys / Innholdsfortegnelse / Spis treści / Índice / Cuprins /
 Содержание / Sadržaj / Obsah / Índice / Innehållsförteckning / İçindekiler

EN	22.5 Fr. Indwelling Voice Prosthesis	6	HU	22,5 Fr. beépített hangprotézisek	79
ZA	22.5 Fr. 留置发音假体	10	IS	22,5 Fr. Inniliggjandi raddventill	84
JA	22.5 Fr. 留置ボイスプロテーゼ	14	IT	Protesi fonatoria fissa da 22,5 Fr.	88
SQ	22,5 Fr. Protezë permanente zëri	18	KO	22.5 Fr. 유치형 인공후두	93
23	جهاز الحنجرة الهوائية المستقرة بمقاس 22.5 على مقياس فرينش	AR	LV	22,5 Fr. Ievadamā balss protēze	97
BG	22,5 Fr. Имплантируема гласова протеза	26	LT	22,5 Fr. Į vidų implantuojamas balso protezas	101
HR	22,5 Fr. Trajna glasovna proteza	32	NO	22,5 Fr. Fast stemmeprotese	106
CS	22,5 Fr. Hlasová protéza	36	PL	22,5 Fr. Stała proteza głosowa	110
DA	22,5 Fr. Indlagt stemmeprotese	41	PT	Prótese vocal fixa de 22,5 Fr.	116
NL	22,5 Ch. Verblijfstemprothese	45	RO	22,5 Fr. Proteză fonatorie permanentă	120
ET	22,5 Fr. Sisemine häälprotees	50	RU	22,5 Fr. Постоянный голосовой протез	125
FI	22,5 Fr. Pysyvä puheproteesi	55	SR	22,5 Fr. Trajna glasovna proteza	131
FR	22,5 Fr. Prothèse vocale à demeure	59	SK	22,5 Fr. Indwelling hlasová protéza	135
DE	22,5 Fr. Verweil-Stimmprothese	64	ES	Prótesis fonatoria fija de 22,5 Fr	140
EL	Μόνιμη φωνητική πρόθεση 22,5 Fr.	70	SV	22,5 Fr. Kvarliggande röstprotes	145
75	22.5 Fr. תותב דיבור פנימי	HE	TR	22,5 Fr. Kalıcı Ses Protezi	149

ENGLISH

INHEALTH 22.5 Fr. INDWELLING VOICE PROSTHESIS

For Secondary placement, Replacement tracheoesophageal puncture (TEP) Procedures

NOTE: This device is not intended to be placed or removed by the patient.

PRODUCT DESCRIPTION

Please refer to the diagrams located at the front of this instruction manual. The voice prosthesis is provided sterile and packaged with one (1) flange introducer and one (1) inserter stick which have all been sterilized by ethylene oxide. Also included in each package are two (2) gel caps. The voice prosthesis is made of medical grade silicone. The flange introducer and inserter stick are both made of thermoplastic (type of plastic) and the gel caps are made of a vegetable-based gelatin. The voice prosthesis contains radiopaque components, which are visible on X-ray. The voice prosthesis consists of a one-way silicone flap valve, an esophageal flange, a body that holds the valve assembly, a tracheal flange, and a safety strap.

INDICATIONS (Reasons to prescribe the device or procedure)

The voice prosthesis is indicated for tracheoesophageal voice restoration following total laryngectomy, when placement, or replacement, of an indwelling prosthesis is performed by a qualified, trained clinician.

CONTRAINDICATIONS (Reasons that make it inadvisable to prescribe the particular device or procedure)

The voice prosthesis is a medical product that should be used only by a qualified clinician with experience and training in its use and care.

WARNINGS AND PRECAUTIONS

Initial sizing and instructions for use must be provided by a qualified clinician trained in the use of this particular device. The clinician is to instruct patients on the use and care of this device and to provide patients with the instructions for use supplied with the product. This voice prosthesis is for single use only. Reuse of this device could result in infection. If the voice prosthesis is not functioning properly, the patient should have it evaluated by a clinician as soon as possible. The voice prosthesis and accessories must be carefully removed from its packaging and handled in a manner to prevent contamination of the device. Handle the device carefully to avoid damage to the device. If there are tears, cracks, or structural damage to the voice prosthesis or its accessories, discontinue use and contact InHealth Technologies Product Complaints. Do not use the device in the event of the sterile packaging being damaged or unintentionally opened prior to use; discard and replace it with a new prosthesis. Use only the appropriately sized gel caps. Do not use petroleum-based products, such as Vaseline® to lubricate the voice prosthesis as these products can damage the device. Changes in the anatomy or medical status of the patient may lead to improper fitting and/or function of the device. Dislodgement or extrusion of the voice prosthesis from the TEP may occur and may result in ingestion, aspiration (inhalation) or tissue damage. A foreign body such as a voice prosthesis in the airway may cause complications such as acute respiratory (breathing) distress and/or respiratory arrest (stop breathing). Always measure the tracheoesophageal puncture (opening made by surgeon between trachea or "windpipe" and esophagus for placement of a voice prosthesis) tract when changing a voice prosthesis to select the proper prosthesis size. A prosthesis, which is too short, may cause tissue necrosis (tissue death) and extrusion (device to protrude out of puncture). Individual patient reactions to the device materials may occur. Consult a clinician immediately if there are indications of tissue edema (swelling) and/or inflammation/infection. Evaluate patients with bleeding disorders or if they are undergoing anticoagulant (keeping blood from clotting) treatment for the risk of hemorrhage (bleeding) prior to placement or replacement of the prosthesis.

VOICE PRODUCTION

To prevent post-operative complications, the patient should not begin speaking with the voice prosthesis until the clinician has indicated that it is safe to do so. The lumen of the voice prosthesis must be kept clear of blockage in order for it to function properly, allowing the patient to voice.

In some patients, the inability to relax the muscles of the throat may account for their inability to speak fluently and to speak with minimal effort. This problem requires professional assessment. Patients requiring postoperative radiation may have transient (brief) interruption of voice in the third or fourth week of treatment. The prosthesis can be left in place or a catheter can be inserted as a stent until radiation treatments are completed.

VOICE PROSTHESIS DISLODGEMENT

Care must be exercised during voice prosthesis insertion or removal to avoid injury to the TEP or accidental displacement of the voice prosthesis, which could result in aspiration (inhalation) of the voice prosthesis into the trachea (windpipe). Should aspiration occur, the patient should attempt to cough the voice prosthesis out of the trachea. Further medical attention may be necessary if coughing the voice prosthesis out is unsuccessful. Confirm gel cap dissolution and deployment of the esophageal flange per procedures provided below to ensure device is securely retained in TEP. Care should be exercised when inserting or removing cleaning devices to avoid accidental displacement of the voice prosthesis, which could result in aspiration of the voice prosthesis. If the voice prosthesis is accidentally dislodged from the puncture, the patient should immediately place a Blom-Singer Puncture Dilator or suitable device of the appropriate diameter, in the puncture to keep it from closing and leaking fluids. The voice prosthesis should be reinserted within 24 hours. Inserting foreign objects other than the Blom-Singer cleaning devices may cause dislodgment of the voice prosthesis or its components and may cause these objects to be swallowed or aspirated.

VOICE PROSTHESIS LEAKAGE

When the valve fails to close completely, a few drops of fluid may pass back through the valve from the esophagus (foodpipe) to the trachea, which may cause coughing. Recurrent leakage should be evaluated by the clinician as leakage could cause aspiration pneumonia. Gentle handling should always be used when cleaning the voice prosthesis to avoid damage, which could cause leakage.

MICROBIAL (MICROSCOPIC ORGANISMS) GROWTH

Most patients have microbial growth in their mouth, throat, or esophagus; this is a commonly occurring condition. Microbial growth deposits on the voice prosthesis may cause valve deformation and failure, i.e., fluid leakage through or around the voice prosthesis and/or an increase in pressure necessary to voice. Replacement of the voice prosthesis may be required.

DEVICE INSERTION OR REMOVAL

The attached safety strap on the voice prosthesis should only be removed after the esophageal flange has been verified to have deployed in the esophagus. Never attempt to insert or reinsert a voice prosthesis that has the safety strap removed. Never remove one voice prosthesis and insert another one without first dilating the TEP and re-measuring the tract to confirm correct voice prosthesis length. The patient should never attempt to insert or remove the voice prosthesis or allow anyone other than a qualified, trained clinician to insert or remove it. The voice prosthesis is not a permanent device and requires replacement periodically. The voice prosthesis may be left in place in the TEP until it ceases to function correctly, that is, until it has persistent leakage, is not providing adequate voice for speech or needs to be resized.

COMPLICATIONS

Although rare, the following complications have been identified to occur with silicone voice prostheses. They include: stoma (opening through neck into trachea or wind pipe) contamination or sepsis, which may require removal of the voice prosthesis and/or appropriate antibiotics; accidental aspiration of the voice prosthesis into the airway, which may require removal by a physician; occasional extrusion of the voice prosthesis, requiring replacement after dilation of the TEP and additional supervision of the stoma care regimen; puncture dilation resulting in leakage of fluids around the voice prosthesis; inflammatory reaction around the puncture site and formation of granulation tissue; dislodgment of the voice prosthesis and subsequent closure of the TE puncture; intractable (uncontrollable) leakage around the voice prosthesis, requiring surgical revision or closure of the puncture; dysphagia (difficulty swallowing), tearing or other damage to the voice prosthesis from improper use; microbial growth deposits causing voice prosthesis leakage or valve incompetence; accidental ingestion of the voice prosthesis into the esophagus.

INSTRUCTIONS FOR USE

NOTE: THIS DEVICE IS NOT INTENDED TO BE PLACED OR REMOVED BY THE PATIENT. THESE INSTRUCTIONS FOR USE ARE FOR THE CLINICIAN ONLY.

The following recommended instructions are provided for primary and secondary/replacement procedures. Please also refer to the diagrams located at the front of this instruction manual.

NOTE: Always measure the TE puncture tract **before** inserting a voice prosthesis to confirm correct voice prosthesis length.

NOTE: Prior to placement of the voice prosthesis, check the valve mechanism to be sure it is intact and working properly. The flap valve should close flat against the seating surface inside the voice prosthesis.

PRIMARY VOICE PROSTHESIS PLACEMENT

The Voice Prosthesis may be placed by the physician at the time of total laryngectomy (primary placement).

SECONDARY VOICE PROSTHESIS PLACEMENT

NOTE: All secondary procedures must only be done using a bright light to illuminate the tracheostoma and puncture.

NOTE: Always dilate TE puncture tract **before** secondary placement of a voice prosthesis. Use of gloves and protective eyewear is recommended.

Prepare the following InHealth Technologies items prior to insertion of the voice prosthesis:

- (1) Voice prosthesis;
- (2) Gel caps;
- (3) Dilator-sizer or sizer with puncture dilator;
- (4) Inserter stick;
- (5) Flange introducer

1. Dilate and Measure (Diagrams 1 & 2)

Please refer to the Blom-Singer Dilation/Sizing System Instructions For Use for complete product details.

2. Voice Prosthesis Preparation

Prepare the voice prosthesis for insertion by utilizing the hand folding method of gel cap insertion. The following steps describe the hand folding method of gel cap insertion.

- 2.1.** Remove the voice prosthesis from the package. Be certain that your hands and the voice prosthesis are completely dry before loading the gel cap. Failure to do so may cause premature dissolution of the gel cap.
- 2.2.** Fold and load the esophageal flange of the voice prosthesis into the gel cap (shorter end of the gel cap) according to the steps below. Discard the gel cap body (longer end of gel cap). (**Diagram 3**)
 - 2.2.1.** Fold the esophageal flange tightly in half toward the center of the device (**Diagram 4**).
 - 2.2.2.** Fold the outer edges of the folded flange tightly against each other (**Diagram 5**) and hold the prosthesis in this folded position.
 - 2.2.3.** Insert the folded portion of the esophageal flange into the gel cap (**Diagram 6**). Gently tuck and push the folded esophageal flange as far as possible into the gel cap.
 - 2.2.4.** Use the curved side of the flange introducer and following the curve of the voice prosthesis body gently tuck the remaining portion of the folded esophageal flange completely into the gel cap (**Diagram 7**). **Do not use any sharp or serrated instruments to avoid damaging the voice prosthesis or gel cap.**
 - 2.2.5.** The esophageal flange should be folded in a forward position and completely enclosed inside the gel cap (**Diagram 8**).

3. Prosthesis Insertion

- 3.1.** Place the voice prosthesis on the inserter stick and lock the safety strap onto the safety peg (**Diagram 9**).
- 3.2.** Caution the patient not to swallow as you remove the dilator-sizer to avoid aspiration of saliva into the trachea prior to inserting the voice prosthesis.
- 3.3.** Apply a light coating of water-soluble lubricant (not petroleum-based) to the tip of the gel-capped end of the voice prosthesis and immediately place the tip of the voice prosthesis in the puncture with the safety strap oriented upwards. Insert the voice prosthesis completely into the dilated TE puncture until the tracheal flange of the voice prosthesis is seated firmly against the posterior tracheal mucosa.

3.4. Hold the voice prosthesis in this position of complete insertion for at least three minutes. This allows ample time for the gel cap to dissolve and release the esophageal flange within the esophagus. Patients should be instructed to swallow their saliva or warm water to facilitate dissolving the gel cap. **Caution:** If the voice prosthesis does not insert easily on the first attempt, do not continue to try to insert it. Instead, re-insert the puncture dilator of the appropriate diameter for a few minutes to dilate the TE tract and then try again.

3.5. Detach the voice prosthesis safety strap from the safety peg of the inserter stick. Place a finger against the safety strap and carefully withdraw the inserter stick from the voice prosthesis with a twisting motion. **Warning:** Caution should be exercised when removing the inserter stick to avoid accidental displacement of the voice prosthesis, which could result in aspiration of the voice prosthesis. Secure the voice prosthesis with medical adhesive tape across the safety strap. It is recommended to orient the prosthesis with the safety strap up.

3.6. Confirm the esophageal flange has been deployed. After insertion of the voice prosthesis in the puncture tract and completion of the 3-minute waiting period for gel cap dissolution and esophageal flange deployment, repeatedly rotate the prosthesis on the inserter stick while it is in the puncture. If the esophageal flange has deployed and has properly seated against the mucosa of the anterior esophageal wall, the voice prosthesis on the inserter stick will rotate easily in the puncture tract. If the gel cap has dissolved and incorrectly deployed the esophageal flange inside the puncture tract rather than inside the lumen of the esophagus, considerable resistance to voice prosthesis rotation will be felt. A second method to confirm gel cap dissolution and esophageal flange deployment is to assess the capacity to produce tracheoesophageal voice. If the gel cap is fully dissolved, the voice prosthesis is correctly positioned and the esophageal flange is deployed, airflow diverted from the trachea by digital occlusion of the tracheostoma will open the flap valve and pass up the esophagus to initiate voice (**Diagram 10**).

3.7. Ensure there is no leakage through the voice prosthesis by having the patient drink water (**Diagram 11**).

4. Safety Strap Detachment

The voice prosthesis is designed to allow the optional detachment of the safety strap following confirmation that the esophageal flange is **securely** positioned within the esophageal lumen against the anterior wall of the esophagus. Removal of the safety strap is an **optional** procedure. Once deployment of the esophageal flange has been confirmed, detach the safety strap from the tracheal flange at the area of reduced strap width, where it meets the tracheal flange, by first locking hemostats onto the strap and then cutting carefully with scissors (**Diagram 12**).

The safety strap should not be detached from the voice prosthesis if:

- Use of a laryngectomy tube or button may potentially dislodge the voice prosthesis;
- A patient has previously been wearing a voice prosthesis that is larger in diameter than the voice prosthesis;
- The diameter of the TE puncture is dilated significantly larger than this voice prosthesis diameter; or,
- The health of the stoma tissue or physical condition of the patient is in question.

Note: Tape the safety strap to the peristomal skin and confirm security of the voice prosthesis retention.

Removal of the Voice Prosthesis

This voice prosthesis is intended to be used for 6 months but this duration should be determined per clinical evaluation and appropriate consultation with a skilled practitioner.

Removal of the voice prosthesis should only be done by grasping the outer tracheal flange of the device securely with a locking hemostat. Pull firmly using caution until the voice prosthesis is fully withdrawn. Insert an appropriately sized dilator-sizer soon after the voice prosthesis is removed and tape it in position for five minutes prior to inserting a new voice prosthesis.

VOICE PROSTHESIS CLEANING AND CARE INSTRUCTIONS

The following information should be made clear to the patient regarding the routine care and cleaning of the indwelling voice prosthesis. The purpose of cleaning the voice prosthesis is to remove blockage of its lumen that impairs the patient's ability to voice while the voice prosthesis is in the patient's TE puncture. Cleaning should be done by using the Blom-Singer® Cleaning Brush

and the Blom-Singer® Flushing Device (sold separately). Please refer to the Cleaning Brush and Flushing Device instructions on how to use the cleaning devices.

Caution: Use **only** the Blom-Singer cleaning devices. **Do not** insert objects other than the Blom-Singer cleaning devices into the voice prosthesis as this may cause damage or dislodgment of the voice prosthesis or its components. Cleaning of the indwelling voice prosthesis should only be done in front of a mirror with a bright light focused directly on the stoma so that the open end of the voice prosthesis is clearly visible. Only use a lint free cloth or tissue to dry the device. Use of non-lint free materials may leave debris, which can be aspirated into the airway. Do not use solvents or petroleum-based products for cleaning or lubricating the device. These materials may damage silicone or cause the device to not work properly.

SPECIAL STORAGE AND/OR HANDLING CONDITIONS

Only the gel caps packaged with the device have any special storage and/or handling conditions. They must be **kept dry** and **not in direct sunlight**.

DISPOSAL INSTRUCTIONS

This product is not biodegradable and must be considered contaminated upon use. Carefully dispose of this device per local guidelines.

ORDERING INFORMATION

USA

Products may be ordered directly from InHealth Technologies. **TELEPHONE:** Toll-Free (800)477-5969 or (805)684-9337, Monday — Friday, 9:30 am — 7:00 pm, Eastern Standard Time. **FAQ:** Toll-Free (888)371-1530 or (805)684-8594. **EMAIL:** order@inhealth.com. **ORDER ON-LINE:** www.inhealth.com. **POST:** InHealth Technologies, 1110 Mark Avenue, Carpinteria, CA 93013-2918, USA, Attention: Customer Service.

INTERNATIONAL

Please contact our customer service department for a distributor referral.

Product Complaints/EU Serious Incidents

If you are dissatisfied with the device or have any questions, please contact productcomplaints@inhealth.com

Telephone: +1-800-477-5969 Fax: +1-888-371-1530

Any serious incident that has occurred in relation to the device should be reported to Freudenberg Medical, LLC as noted above and the competent authority of the EU Member State in which the user and/or patient is established.

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SYMBOLS GLOSSARY (USA only)

SYMBOL	STANDARD No.	SYMBOL TITLE & REF No.	MEANING OF SYMBOL
	EN ISO 15223-1:2016	Manufacturer; 5.1.1	Indicates the medical device manufacturer, as defined in EU Directives or Regulations
	93/42/EEC EU 2017/745	CE Marking; 93/42/EEC Annex XII, EU 2017/745 Annex V	CE marking of conformity
	EN ISO 15223-1:2016	Authorized representative in the European Community; 5.1.2	Indicates the authorized representative in the European Community
	EN ISO 15223-1:2016	Catalogue number; 5.1.6	Indicates the manufacturer's catalogue number so that the medical device can be identified
	EN ISO 15223-1:2016	Batch code; 5.1.5	Indicates the manufacturer's batch code so that the batch or lot can be identified
	EN ISO 15223-1:2016	Use-by date; 5.1.4	Indicates the date after which the medical device is not to be used
	EN ISO 15223-1:2016	Date of manufacture; 5.1.3	Indicates the date when the medical device was manufactured
	EN ISO 15223-1:2016	Caution; 5.4.4	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself

SYMBOLS GLOSSARY (USA only)

SYMBOL	STANDARD No.	SYMBOL TITLE & REF No.	MEANING OF SYMBOL
	EN ISO 15223-1:2016	Do not use if package is damaged; 5.2.8	Indicates a medical device that should not be used if the packaging has been damaged or opened
	EN ISO 15223-1:2016	Consult instructions for use; 5.4.3	Indicates the need for the user to consult the instructions for use
	EN ISO 15223-1:2016	Do not re-use; 5.4.2	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure
	EN ISO 15223-1:2016	Non-sterile; 5.2.7	Indicates a medical device that has not been subjected to a sterilization process
	EN ISO 15223-1:2016	Sterilized using ethylene oxide; 5.2.3	Indicates a medical device that has been sterilized using ethylene oxide
	EN ISO 15223-1:2016	Keep away from sunlight; 5.3.2	Indicates a medical device that needs protection from light sources
	EN ISO 15223-1:2016	Keep dry; 5.3.4	Indicates a medical device that needs to be protected from moisture
	ISO 7000	Packaging unit; 2794	To indicate the number of pieces in the package

INHEALTH

TECHNOLOGIES®

UNITED STATES

1110 Mark Avenue
Carpinteria, CA 93013-2918 USA
Toll-Free: 800.477.5969
Toll-Free Fax: 888.371.1530

INTERNATIONAL

Tel: +1.805.684.9337
Fax: +1.805.684.8594
E-mail: order@inhealth.com

ORDER ONLINE

www.inhealth.com

Freudenberg Medical, LLC
1110 Mark Avenue
Carpinteria, CA 93013-2918 USA
www.inhealth.com

CE 0344

EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands