












LARYVOX® UNIEK

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
PICTROGRAM LEGEND

Pictogrammes listed below you'll find on the product packaging, if applicable.

en

	Catalogue number
	Batch code
	Use by
	Date of manufacture
	Content (in pieces)
	Consult instructions for use
	Keep away from sunlight
	Keep dry
	Do not reuse
	Medical devices
	Sterilized using ethylene oxid
	Do not resterilize
	Do not use if package is damaged
	Simple sterile barrier system
	CE marking identification number of the notified body
	Manufacturer

LARYVOX® UNIEK

1. FOREWORD

These instructions are for **LARYVOX® UNIEK**.

These instructions for use serve to inform the physician, nursing staff and patient/user in order to ensure technically correct handling.

The selection, use and insertion of the products must be carried out by a trained physician.

Please read the instructions for use carefully before using the product for the first time!

Keep them in an easily accessible place for future reference.

Please keep the package for as long as you use the product.

It contains important information on the product!

2. SAFETY INFORMATION

Both the intraoperative placement as well as the routine replacement of the voice prosthesis may only be performed by qualified medical specialists trained in accordance with local and national guidelines.

The user must first have been trained by qualified medical personnel in the safe handling and application of the **LARYVOX® UNIEK**.

LARYVOX® UNIEK is a single-patient product and thus intended exclusively for use in a single patient.

Disinfection, re-sterilisation or reuse can impair the safety and function of the product and are therefore not permissible!

Reuse with other patients and preparation for re-use are not permissible.

The time of the first attempt to speak after the operation (initial installation) falls into the physician's sphere of responsibility in order to avoid possible complications.

Do not perform any type of adjustment, repair or alteration on the voice prosthesis.

In the event of damages, contact the manufacturer.

In general, avoid exerting external pressure on the voice prosthesis using accessories, e.g. tracheal cannulas or tracheostomy buttons, in order to prevent dislocation or the development of tissue damage.

Before replacing the voice prosthesis, take into account the patient's blood coagulation values.

If the patient takes anticoagulants or is suffering from a blood coagulation disorder, there is an increased risk of haemorrhage.

The voice prosthesis must be cleaned at regular intervals (see Chapter 9. HYGIENE INSTRUCTIONS) in order to preserve the functionality and voice function for the patient.

Incorrect handling of accessories can lead to damage of the voice prosthesis.

3. INTENDED PURPOSE

The **LARYVOX® UNIEK** is a long-term voice prosthesis made of medical silicone for voice rehabilitation after a laryngectomy.

The **LARYVOX® UNIEK** is intended for long-term application.

The voice prosthesis prevents closure of the tracheoesophageal shunt canal. It enables air to flow into the pharynx and thus allows speech.

The voice prosthesis seals the tracheoesophageal shunt canal during swallowing.

4. INDICATION

The use of the **LARYVOX® UNIEK** voice prosthesis is indicated for laryngectomised patients with a tracheoesophageal shunt canal.

It can be introduced into the tracheo-oesophageal shunt channel/fistula either directly in the scope of tumour surgery (primary puncture), or at a later time (secondary puncture).

In the event of loss of function during use, the voice prosthesis must be replaced promptly.

The **LARYVOX® INSERTER** facilitates the anterograde placement/change of a voice prosthesis.

A voice prosthesis cannot remain permanently in the shunt channel/fistula, but must be replaced at regular intervals, depending on how it is functioning.

5. CONTRAINDICATION

Do not use if the patient is allergic to the material.

6. COMPLICATIONS AND ADVERSE REACTIONS

Please contact your treating physician if any of the following occur:

- persistent cough and/or abnormal discolouration of the secretion (mainly in connection with food and fluid intake)
- alteration / deterioration of the voice quality
- pain, swellings and traces of blood after cleaning

We recommend routine examinations on a regular basis in consultation with your doctor.

General changes of the state of your tissue or of the general health condition can impair the seating and/or the function of the voice prosthesis.

It is recommended to have the voice prosthesis and the shunt channel/fistula checked on a regular basis by a medical professional.

If cleaning accessories are used, then the respective instructions for use must be observed.

Dislocation of the voice prosthesis

If the voice prosthesis is no longer properly in place in the shunt channel/fistula and protrudes markedly from the shunt channel/fistula into the trachea, or the tracheal flange is only barely still visible, then a dangerous change of position has taken place, which must be checked/corrected immediately by a physician.

The oral intake of fluids and/or nutrition must be absolutely avoided.

Swallowing of the voice prosthesis

If the voice prosthesis is no longer located in the shunt channel/fistula, then the voice prosthesis may possibly have been swallowed or aspirated.

A physician must be consulted immediately.

The oral intake of fluids and/or nutrition must be absolutely avoided.

Ingestion of the voice prosthesis into the stomach may be possible.

Foreign matter which accidentally enters the stomach normally passes through the digestive tract and is then excreted.

If severe symptoms occur, such as haemorrhages, ileus and/or sepsis, then it must be assumed that the voice prosthesis was not excreted.

A surgical intervention to remove the voice prosthesis must be considered.

If the voice prosthesis that has been swallowed is situated in the lower section of the oesophagus, then the voice prosthesis can be located by means of an X-ray and removed by oesophagoscopy.

Aspiration of the voice prosthesis

If the voice prosthesis is no longer located in the shunt channel/fistula, it is possible that the voice prosthesis may have been inhaled.

Symptoms of possible aspiration are:

- Foreign matter sensation in the airways with persistent urge to cough
- Breath sounds
- Difficulty breathing
- Urge to gag

A physician must be consulted immediately for removal of the aspirated voice prosthesis.

The oral intake of fluids and/or nutrition must be absolutely avoided.

If the condition is not treated or treatment is delayed, further complications (pneumonia, atelectasis, bronchitis, lung abscess, bronchopulmonary puncture and asthma attack) may arise.

Measures for avoiding dislocation / aspiration / swallowing

A daily inspection and assessment of the voice prosthesis with regard to its position in the shunt channel/fistula as well as alterations of the shunt channel/fistula is absolutely required in order to detect and avoid dislocations.

Therapeutic appliances such as for instance tracheostomy tubes or tracheostomy buttons must be selected individually and must not exert any pressure on the voice prosthesis and the shunt channel/fistula.

Leakage of the voice prosthesis

The leakage of fluids with intrusion of saliva or food debris from the oesophagus into the trachea, accompanied by the urge to cough, is a direct indication for an inspection of the voice prosthesis.

By the patient: check during fluid intake, whether the leakage occurs centrally through the voice prosthesis (transprosthethically) or at the side out of the shunt channel/fistula (periprosthethically).

Procedure:

1. in the event of transprosthethic leakage

- Clean the voice prosthesis with the **LARYVOX® DUO BRUSH** in order to remove possible food debris / incrustations.
- Check the leakproofness of the voice prosthesis by means of a drinking test.
- If the leakage persists, insert the **LARYVOX® PLUG** and promptly consult your treating physician.

2. periprosthethic leakage

- Consult your doctor immediately.
Avoid any further food / fluid intake.
- The voice prosthesis may have to be replaced.

By the physician: visual inspection of the voice prosthesis while the patient performs a drinking test and, if necessary, replacement of the voice prosthesis used.

In the event of chronic periprosthethic escape of fluids, surgical measures or closure of the shunt channel/fistula may become necessary.

Voice / swallowing impairments

Check the shunt valve for patency (e.g. using the **LARYVOX® DUO BRUSH** cleaning brush).

If voice / swallowing impairments persist, consult your treating physician in order for instance to have the voice prosthesis length checked.

Candida infection

Candida (yeast fungus) infection is a common phenomenon that influences the function of the voice prosthesis.

Radiotherapy, gastro-oesophageal reflux and eating habits influence the timeline of colonisation of the voice prosthesis with candida and the resulting transprosthethic defect.

Minor haemorrhages

Haemorrhages on the shunt channel/fistula are not uncommon when performing routine replacement of the voice prosthesis.

Intensified haemorrhage may occur under administration of anticoagulants.

Granulation tissue

Granulation tissue may occasionally form around the shunt channel/fistula.

If this is the case, removal by electrical or chemical means or by means of a laser cutter should be considered.

Necrosis formation

Formation of necroses on the mucosa on the tracheal or oesophageal side may occur due to too short voice prosthesis lengths.

Scar formation

Excessive scar formation is promoted by too short voice prosthesis.

The question, whether the scar formation can be removed or whether a longer voice prosthesis must be used, must be weighed carefully.

Further complications

Infections, contaminations, mucosal oedemas, proliferations/granulations, change of shunt channel/fistula length after exposure to radiation, inflammations and scar tissue at the puncture site, enlargement of the diameter of the shunt channel/fistula.

If complications occur, please contact your treating physician immediately.

7. PRODUCT DESCRIPTION

The package contains 1 **LARYVOX® UNIEK**, packed sterile and sterilised with ethylene oxide (EO), 1 **LARYVOX® INSERTER UNIEK** placement instrument and 1 **LARYVOX® DUO BRUSH** cleaning brush.

The **LARYVOX® UNIEK** voice prosthesis is used for voice rehabilitation following removal of the larynx (laryngectomy).

This prosthesis is a placeholder for a surgically created channel between the trachea and the oesophagus (tracheo-oesophageal shunt channel/fistula) that serves for voice rehabilitation.

The **LARYVOX® UNIEK** is a product made of medical-grade silicone.

The design comprises a tracheal-side flange (tracheal flange) and an oesophageal-side flange (oesophageal flange) with a valve flap on the oesophageal side and a safety thread on the tracheal flange for placement in the shunt channel/fistula (Figure 1).

en

If the tracheostoma is closed during exhalation, the airflow via the **LARYVOX® UNIEK** is directed towards the oesophagus, causing vibrations of the different tissue structures to create your individual voice.

The intrusion of food or fluids into the windpipe (trachea) is prevented by the valve flap.

The **LARYVOX® UNIEK** is available in various shaft lengths: 5, 6, 7, 8, 9, 10, 11, 13.

The shaft diameter is 7mm (21 Fr.) or 8mm (24 Fr.), in order to enable precise adjustment to different shunt channel/fistula sizes.

The **LARYVOX®** voice prosthesis can remain in the shunt channel/fistula during a magnetic resonance tomography (MRT).

Radiation doses of more than 70 Gy within a time period of 7 weeks increase the risk of tissue necroses in and around the region of the shunt channel/fistula.

The creation of a shunt channel/fistula for voice rehabilitation by means of voice prosthesis is not advisable under these circumstances.

The **LARYVOX® UNIEK** is made of a radiocontrasting material and can be depicted in radiodiagnostics.

A loss of voice may occur during the course with radiotherapy.

This is not the result of a malfunction of the voice prosthesis, but rather a side effect of the associated tissue irritation during radiotherapy.

If the voice prosthesis continues to be seated correctly, then the voice prosthesis can remain in place.

ACCESSORIES:

LARYVOX® INSERTER (REF 25925)

The **LARYVOX® INSERTER UNIEK** is used for the anterograde placement/replacement of a **UNIEK** voice prosthesis.

LARYVOX® DUO BRUSH (REF 29110)

Cleaning brush for daily cleaning of the inner channel of the voice prosthesis.

LARYVOX® PLUG (REF 25800)

The **LARYVOX® PLUG** sealing plug serves for short-term use with a voice prosthesis in order to prevent transprosthetic fluid leakage (centrally through the voice prosthesis).

8. APPLICATION

WARNING!

The voice prosthesis may only be changed by qualified medical specialists who have been trained in accordance with local and national guidelines.

WARNING!

Never use the LARYVOX® UNIEK or the LARYVOX® INSERTER if these items show visible faults or defects.

Contact the manufacturer in such cases.

WARNING!

A new voice prosthesis must only be inserted after the length of the shunt channel/fistula has been checked and the correct length of prosthesis has been selected.

It is recommended to use a special measuring instrument to check the length of the shunt channel/fistula in order to determine the correct shunt channel / fistula length prior to placement.

CAUTION!

The voice prosthesis must always be placed with care in the shunt channel/fistula.

Pay attention to the correct position of the voice prosthesis.

The **LARYVOX® INSERTER** is used for anterograde placement / replacement of a voice prosthesis.

The voice prosthesis is firmly attached to the placement instrument with the **LARYVOX® INSERTER** and is then introduced into the existing shunt channel/fistula using the appropriate applicator.

Please also observe the notes / warnings in the instructions for use for the **LARYVOX® INSERTER**.

It is advisable to always perform the replacement of the voice prosthesis with optimal illumination of the tracheostoma and the shunt channel/fistula.

The placement of the voice prosthesis must without fail be carried out under hygienic conditions (gloves, eye protection etc.).

Before removing the voice prosthesis to be replaced, check its seating in the shunt channel/fistula.

The correct length of the voice prosthesis is of major importance for avoiding complications such as periprosthetic leakage, pressure necroses or the like.

Replacing the voice prosthesis

Carefully check the sterile packaging to ensure that it has not been tampered with or damaged.

Do not use the product if the packaging has been damaged.

Observe the use-by/expiry date.

Do not use the product after this date.

PLEASE NOTE!

Based on experience, 2 - 3 placement attempts with the **LARYVOX® INSERTER** are possible in the scope of one application.

CAUTION!

Use only and exclusively original **LARYVOX®** accessories.

PLEASE NOTE!

Under no circumstances should silicone spray or a fatty cream be used as an additional lubricant.

PLEASE NOTE!

Check the integrity of the valve flap before replacement.

1. Fit the voice prosthesis with the tracheal flange onto the head of the placement instrument (Figure 2).
2. Fix the safety thread of the voice prosthesis in place on the placement instrument at position A (Figure 3).
3. Check the secure seating of the voice prosthesis on the placement instrument.
4. To achieve better lubrication of the applicator for advancing the voice prosthesis, it is recommended to use the **FAHL® OPTIFLUID® LUBRICANT GEL** included in delivery.
5. Use the placement instrument to push the voice prosthesis forward by approx. 1 cm into the slotted opening of the applicator (Figure 4).
Marker positions 1 and 2 are shown on the placement instrument.
6. Now guide the placement instrument with the voice prosthesis through the applicator up to the line of marker position 1 (Figure 5).
7. Cautiously introduce the tip of the applicator into the shunt channel (Figure 6).
8. Now slowly advance the voice prosthesis through the applicator up to the line of marker position 2 (Figure 7).
The oesophageal flange opens in the lumen of the oesophagus.
9. Now retract the applicator carefully via the placement instrument (Figure 8).
The tracheal flange opens in the lumen of the trachea.
If the voice prosthesis cannot be inserted at the first attempt, the placement procedure should be repeated after a brief pause.
10. The correct position of the voice prosthesis is verified by rotating the voice prosthesis with the placement instrument or by pulling slightly (Figure 9).
11. Rotate the voice prosthesis into position such that the safety thread is in the 12 o'clock position.
12. Fix the safety thread in place with an atraumatic clamp and remove the placement instrument (Figure 10).
Never grasp the voice prosthesis with the instrument, as this can damage or dislocate the voice prosthesis.
Apply only light resistance when removing the placement instrument.

WARNING!

Proceed with particular care when doing this in order to avoid dislocation, even up to aspiration of the voice prosthesis.

13. Perform a phonation test.
14. Check the leakproofness of the voice prosthesis with a swallowing test.
15. If the voice prosthesis is in the correct position, cut the safety thread flush with the tracheal flange (Figure 11).

WARNING!

The safety thread on the tracheal flange must only be severed after correct positioning of the voice prosthesis in the shunt channel/fistula.

Placement of the voice prosthesis without safety thread is prohibited.

Removing the voice prosthesis

WARNING!

Fixation of the tracheal flange with a fixable atraumatic clamp is absolutely required for the safe removal of the voice prosthesis.

The voice prosthesis can be withdrawn completely from the tracheo-oesophageal shunt channel by carefully pulling on the tracheal flange.

The currently used voice prosthesis can be used as orientation for determining the correct shunt valve length.

1 - 2 mm clearance between the mucosa and the tracheal flange constitute an optimal fit.

Smaller or greater distances require an adjusted size selection of the new voice prosthesis.

9. HYGIENE INSTRUCTIONS

Clean the voice prosthesis at regular intervals in order to preserve the functionality and voice function.

CAUTION!

Cleaning activities of any kind on and in the voice prosthesis must be carried out with particular diligence and caution, exerting as little pressure as possible, in order to avoid dislocation / aspiration / swallowing of the voice prosthesis.

CAUTION!

Before using the cleaning accessories, such as for instance **LARYVOX® DUO BRUSH**, or before inserting the **LARYVOX® PLUG**, hands must be cleaned with water and soap in order to avoid contaminations and infections.

Use only and exclusively the **LARYVOX® DUO BRUSH** cleaning brush for daily cleaning of the voice prosthesis.

CAUTION!

The voice prosthesis must be cleaned only and exclusively within the shunt channel/fistula.

Never attempt to remove the voice prosthesis from the shunt channel/fistula for cleaning.

Avoid inserting the cleaning brush incorrectly with slanting/tilting of the voice prosthesis, which could lead to dislocation of the voice prosthesis.

Do not use the **LARYVOX® DUO BRUSH** cleaning brush if the shunt channel/fistula is greatly widened, such that the secure hold of the voice prosthesis in the tissue is not ensured.

Please also observe the notes / warnings in the instructions for use for the **LARYVOX® DUO BRUSH**.

Check the **LARYVOX® DUO BRUSH** cleaning brush for possible damage before each use.

Only undamaged cleaning brushes must be used for cleaning the voice prosthesis.

Do not perform any type of adjustment, repair or alteration on the brush or compatible accessories.

If products are damaged, they must be disposed of properly.

Cleaning the voice prosthesis

Moisten the brush/sponge with clean water before use.

Perform the cleaning in front of a mirror under optimal lighting conditions so that the voice prosthesis is clearly visible (Picture 12).

If necessary, use the **LARYVOX® CLIP** with a tracheostoma lamp.

A stoma lamp can be affixed to the cleaning brush by means of the **LARYVOX® CLIP** in order to provide good illumination of the voice prosthesis for cleaning (Figure 13).

Gently and cautiously introduce the cleaning brush into the inner channel of the voice prosthesis (Figure 14).

Clean the voice prosthesis by gentle forward-and-back movements accompanied by simultaneous gentle turning movements (Figure 15).

To optimise the accessibility of the voice prosthesis, the cleaning brush can be bent individually at the flexible brush body (not at the brush wire!) (Figure 16).

Resterilisation of the voice prosthesis is impermissible.

10. STORAGE

This product should be stored in a dry environment away from sunlight and/or heat.

11. SERVICE LIFE

A voice prosthesis cannot remain permanently in the shunt channel/fistula, but needs to be replaced at regular intervals, depending on its proper functioning.

The service life of voice prostheses generally depends on the individual circumstances and the patient's metabolic situation.

Accordingly, it is not possible to make a clear prognosis of the service life.

Particularly the biofilm colonisation of the voice prostheses with bacteria and yeast fungi have an impact on the service life of the product.

Candida infection is a frequently occurring phenomenon which influences the function of the voice prosthesis.

Radiotherapy, gastro-oesophageal reflux and eating habits influence the timeline of colonisation of the voice prosthesis with candida and the resulting trans-prosthetic defect.

Leakages constitute the main indications for required replacement.

Leakages can occur through the voice prosthesis, or around the voice prosthesis.

Anatomic changes of the shunt channel/fistula, infections or other medical indications, too, can necessitate the replacement of the voice prosthesis.

Use beyond 12 months is at the discretion of the treating physician.

The **LARYVOX® UNIEK VOICE PROSTHESIS** is intended for one-time use.

12. DISPOSAL

The product must only be disposed of in accordance with the national regulations.

13. LEGAL NOTICES

The manufacturer Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for functional deficiencies, injuries, infections, and/or other complications or adverse events caused by unauthorised product alterations or improper use, care, and/or handling.

In particular, Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for damages caused by modifications to the product or by repairs, if these modifications or repairs were not carried out by the manufacturer himself.

This applies – to the extent permitted by law – both to damages to the products themselves caused thereby as well as to any consequential damages caused thereby.

When using the product beyond the period of use specified under Section 11, and/or use, treatment, maintenance (cleaning, disinfection) or storage of the product in non-compliance with the specifications laid down in these instructions for use, Andreas Fahl Medizintechnik-Vertrieb GmbH will, as far as legally permissible, be free of any liability, including liability for defects.

Should a serious incident occur in connection with this product of Andreas Fahl Medizintechnik-Vertrieb GmbH, then this must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient are domiciled.

Sale and delivery of all Andreas Fahl Medizintechnik-Vertrieb GmbH products is carried out exclusively in accordance with our General Terms and Conditions of Business (GTC) which can be obtained directly from Andreas Fahl Medizintechnik-Vertrieb GmbH.




The manufacturer reserves the right to make changes to the product at any time without notice.

LARYVOX® is a trademark and brand of Andreas Fahl Medizintechnik-Vertrieb GmbH, Cologne registered in Germany and the EU member states.

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LARYVOX® är ett registrerat varumärke i Tyskland och EU som tillhör Andreas Fahl Medizintechnik-Vertrieb GmbH, Köln.