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
# PICTROGRAM LEGEND

en

Pictogrammes listed below you'll find on the product packaging, if applicable.

Catalogue number

Batch code

Use by

Date of manufacture

Content (in pieces)

Consult instructions for use

Keep away from sunlight

Keep dry

Do not reuse

Medical devices

Sterilized using ethylene oxid

Do not reesterilize

Non-sterile

Do not use if package is damaged

Single sterile barrier system

CE marking

CE marking identification number of the notified body

Manufacturer

# HME / ARTIFICIAL NOSE

## 1. FOREWORD

These instructions for use are valid for Artificial Noses/HME (Heat-Moisture Exchanger) (s. cover sheet).

These instructions for use serve to inform the physician, nursing staff and patient/user in order to ensure technically correct handling.

When used for the first time, the products must be selected, used and inserted by appropriately trained physicians or medical professionals/medical device consultants.

Please read the instructions for use carefully before using the product!

Keep them in an easily accessible place for future reference.

Please keep the package for as long as you use the product.

## 2. SAFETY INFORMATION

All HMEs are intended for single patient use only and must not be cleaned or disinfected.

HMEs must under no circumstances be inserted directly into the tracheostoma.

If the HME is impermeable (due to secretion) or shows signs of impermeability, it must be removed and disposed of immediately.

Care must be taken to ensure that the airflow is not blocked by clothing or the like (risk of respiratory distress) to enable unhindered breathing.

Do not aspirate secretion through the valve opening under any circumstances.

If the filter material has slipped in the housing, the product must be replaced and disposed of.

The HMEs must not be fixed to the tracheostomy tube while using a nebuliser or an atomiser.

Do not fill any water into the HME.

When used in combination with tracheostomy tubes and base plates, compatibility and a secure fit must be checked with utmost care.

HMEs must not be disassembled and then reassembled as this will impair their function.

Do not use on patients with too low a respiratory volume!

Respiratory resistance and additional dead space can lead to a build-up of carbon dioxide and cause acute respiratory distress.

Do not grease connections, tubes or other accessories in order to avoid fire hazard.

Smoking during oxygen therapy is dangerous and can lead to severe fire-induced injury.

Only use lotions or ointments labelled as oxygen-compatible to avoid the risk of fire or burns.

There is a fire hazard in connection with oxygen devices and oxygen therapies.

Do not use in the vicinity of sparks or open flames.

Patients who have not previously used HMEs may experience increased respiratory resistance during initial use.

It is recommended to increase the application period slowly.

## 3. INTENDED PURPOSE

HMEs serve to humidify the airways and warm the air one breathes.

They reduce the formation of viscous secretions in the lungs and prevent the trachea from drying out.

HMEs filter the breathing air, thereby reducing intrusion of particles into the patient's airways.

The use of an HME creates the respiratory resistance that is necessary to maintain physiological lung function.

Variants with an O<sub>2</sub> connection also allow for the supply of additional oxygen.

## 4. INDICATION

en

HMEs are suitable for spontaneously breathing patients with a tracheostoma.

The products **HUMIDOSTOM® JUNIOR (REF 46850)** and **HUMIDOSTOM® MINI (REF 46870)** are intended for use in newborns and children.

## 5. CONTRAINDICATION

The HME should be removed if respiratory resistance is too high.

The use in patients with severe exsiccosis (dehydration), excessive secretion and water accumulation in the lungs and airways, as well as very viscous secretion and haemoptysis is not permitted.

This product is not suitable for patients with low respiratory volumes as the additional dead space may lead to CO<sub>2</sub> retention.

It must not be used outside the recommended tidal volume range, as increased dead space can lead to CO<sub>2</sub> retention if the tidal volume is too low and to insufficient humidification if the tidal volume is too high.

For persons with impaired consciousness, patients with limited motor function and/or persons requiring a high level of care, and who are not in a position to remove the product by themselves, monitoring by trained specialist personnel is required to detect excessive breathing resistance in time and to initiate the appropriate measures to safeguard the airways.

## 6. COMPLICATIONS AND ADVERSE REACTIONS

Physical exertion can lead to an increased effort to breathe and shortness of breath due to respiratory resistance.

In this case, it is recommended to withdraw the HME for a short period of time and to train using the HME long-term to achieve better acceptance in the future.

## 7. PRODUCT DESCRIPTION

HMEs consist of a foam filter (**REF 46440** consists of a paper filter) and a plastic housing.

The plastic housing is equipped on the patient side with a central opening with 15 mm internal diameter, thereby allowing for connection to a 15 mm connector.

**REF 46860**: Connection to **TRACHLINE®** tracheostomy tube via bayonet connector.

**REF 46150, 46155, 46160**: The housing cover has a grid design to ensure an unobstructed air flow during coughing.

**REF 46150, 46155, 46450, 46456, 46457, 46470**: an oxygen connecting piece fitted additionally to the side of the housing enables attaching a connecting tube, thereby allowing combination with an oxygen source.

Here, a distinction is made between a 5 mm plug-in connection (**REF 46150, 46450, 46456, 46470**) and a 6 mm O<sub>2</sub> ISO connection (**REF 46155, 46457**).

**REF 46440**: The filter medium in the transparent plastic housing consists of hydrophilic, moisture-binding cellulose.

**REF 46440, 46470**: The product was packaged singly and sterilised with ethylene oxide (EO). In addition, an adapter with oxygen connection can be ordered as a separate accessory (**REF 46441**).

**REF 46460**: The external diameter of the connector is 22 mm, and is thus compatible with all commonly used 22 mm receptacle systems.

The functional multi-adapter features both a central opening for attachment onto a tracheostomy tube with 15 mm connector, as well as a 22 mm combi-adapter.

**REF 46450, 46456, 46457, 46460, 46470**: The overpressure valve in the lid of the housing prevents an excessive pressure build-up and allows secretions to be coughed up without having to remove the HME.

**REF 46850:** HME with low weight, specifically for children and newborns.

**REF 46860:** The plastic housing is equipped on the patient side with a bayonet lock and thereby ensures connection to a **TRACHLINE®** tracheostomy tube.

In addition, it features an oxygen connection (ISO connection, 6 mm).

**REF 46870:** HME for children with side-mounted oxygen connection (ISO connection, 6 mm).

	Oxygen connector	Overpressure valve	Grid-shaped housing cover	Lateral grid-shaped housing cover	sterile
<b>REF 46150</b>	x		x		
<b>REF 46155</b>	x		x		
<b>REF 46160</b>			x		
<b>REF 46440</b>				x	x
<b>REF 46450</b>	x	x			
<b>REF 46456</b>	x	x			
<b>REF 46457</b>	x	x			
<b>REF 46460</b>		x			
<b>REF 46470</b>	x	x		x	x
<b>REF 46850</b>					
<b>REF 46860</b>	x				
<b>REF 46870</b>	x				

## 8. APPLICATION

1. Check whether the packaging and the product are in a faultless condition.
2. Check the use-by/expiry date.  
Do not use the product after this date.
3. Check the stated tidal volume.
4. Remove the existing HME by holding the tracheostomy tube/base plate in place with two fingers and pulling the HME off gently (Fig. 3).
5. Mount the new HME with the central 15 mm opening onto the 15 mm connector of the tracheostomy tube/endotracheal tube or the 22 mm connector of the base plate by applying slight pressure (Fig. 1).
6. For an HME with 22 mm combination adapter (**REF 46460**):  
mount the HME on the matching connector of the tracheostomy tube/base plate by applying slight pressure, the housing will snap noticeably into place through the outer mounting rim (Fig. 4).  
**REF 46860:** Mount the HME onto the 15 mm connector of the tracheostomy tube with the bayonet lock.  
Fixate the HME in the tracheostomy tube with a slight anticlockwise rotation (on the patient's side when inserted), thus closing the bayonet lock.  
**REF 46150, 46450, 46455, 46450, 46456, 46457, 46470:** For additional administration of oxygen, a connecting tube can be connected to the O<sub>2</sub> connection.  
To this end, push the connecting tube onto the O<sub>2</sub> connection.

## 9. HYGIENE INSTRUCTIONS

HMEs are single-patient products and intended for single use only.

They must not be cleaned or disinfected, as use of disinfectant can render the HME unusable.

In addition, bacterial or fungicidal colonisation of the filter material can lead to an increased risk of infection.

Cleaning, disinfection or (re-)sterilisation as well as use longer than 24 hours can impair the safety and function of the product and are therefore not permissible!

## 10. STORAGE

HMEs should be stored in a dry environment away from sunlight and/or heat.

## 11. SERVICE LIFE

The maximum period of use / wearing period is 24 hours (Fig. 2).

HMEs are even to be changed several times a day if required, for instance if the respiratory resistance has increased to an unacceptable level due to accumulation of secretions.

## 12. DISPOSAL

The product must only be disposed of in accordance with the national regulations.

## 13. LEGAL NOTICES

The manufacturer Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for functional deficiencies, injuries, infections, and/or other complications or adverse events caused by unauthorised product alterations or improper use, care, and/or handling.

In particular, Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for damages caused by modifications to the product or by repairs, if these modifications or repairs were not carried out by the manufacturer himself.

This applies both to damage caused to the products themselves as a result thereof and to all consequential damage caused thereby.

When using the product beyond the period of use specified under Section 11, and/or use, treatment, maintenance (cleaning, disinfection) or storage of the product in non-compliance with the specifications laid down in these instructions for use, Andreas Fahl Medizintechnik-Vertrieb GmbH will, as far as legally permissible, be free of any liability, including liability for defects.

Should a serious incident occur in connection with this product of Andreas Fahl Medizintechnik-Vertrieb GmbH, then this must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient are domiciled.

Sale and delivery of all Andreas Fahl Medizintechnik-Vertrieb GmbH products is carried out exclusively in accordance with our General Terms and Conditions of Business (GTC) which can be obtained directly from Andreas Fahl Medizintechnik-Vertrieb GmbH.

The manufacturer reserves the right to make changes to the product at any time without notice.

**HUMIDOSTOM®**, **HUMIDOTRACH®**, **HUMIDOFIX®**, **HUMIDOTWIN®**, **HUMIDUAL®** and **TRACH-LINE®** are registered trademarks of Andreas Fahl Medizintechnik-Vertrieb GmbH, Cologne in Germany and the EU member states.

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