





MUCOPROTECT®

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PICTOGRAM LEGEND

Pictogrammes listed below you'll find on the product packaging, if applicable.

Catalogue number

Batch code

Use by

Consult instructions for use

Keep away from sunlight

Keep dry

Medical devices

CE marking identification number of the notified body

Manufacturer

en

MUCOPROTECT® CUFF PRESSURE GAUGE

1. FOREWORD

This instruction applies for the **MUCOPROTECT®** cuff pressure gauge.

These instructions for use serve to inform the physician, nursing staff and patient/user in order to ensure technically correct handling.

When used for the first time, the products must be selected, used and inserted by appropriately trained physicians or medical professionals/medical device consultants.

Please read the instructions for use carefully before using the product for the first time!

Keep them in an easily accessible place for future reference.

Please keep the package for as long as you use the product.

It contains important information on the product!

2. SAFETY INFORMATION

When using the cuff pressure gauge for monitoring, the cuff pressure must be checked and adjusted by the user every 15 minutes.

Upon attachment of the connecting tube to the control balloon of the tracheotomy tube/tracheal tube, the connecting tube opens the balloon valve.

When removing the cuff pressure gauge, the cuff pressure gauge must be removed together with the connecting tube, as otherwise the air will escape uncontrollably from the cuff of the tracheotomy tube/tracheal tube.

Patients must be trained by medical professionals on how to use safely.

The cuff pressure gauge is not compatible with MRT!

The connecting tube is a single-patient product and thus intended exclusively for use in a single patient.

Re-use of any kind in other patients and thus also reconditioning for re-use in other patients is not allowed!

3. INTENDED PURPOSE

The cuff pressure gauge is used for ventilation as well as for pressure regulation and control of a tracheostomy tube/ tracheal tube with cuff.

4. INDICATION

Tracheostomy tubes with cuff are always indicated after tracheostomy has been performed for any cause if the cannula needs to be sealed against the wall of the trachea.

The cuff itself is inflated via a tube with one-way valve and pilot balloon.

The cannula and the cuff must be tested for leakages directly before and after every insertion and thereafter at regular intervals.

5. COMPLICATIONS AND ADVERSE REACTIONS

It is recommended that the balloon internal pressure be kept above 20 cm H₂O because of the risk of aspiration or pneumonia, but also below the value of 30 cm H₂O because of the risk of ischemia of the tracheal mucosa.

The recommended values may be varied by doctor's orders.

6. PRODUCT DESCRIPTION

Technical data

Measuring range:	0 – 100 cm H ₂ O
Storage conditions:	-30°C to 70°C, 15 – 85 % relative humidity
Operating conditions:	10°C – 40°C, 15 – 85 % relative humidity
The measurement accuracy is ± 2 cm H ₂ O for the entire display range.	1 cm H ₂ O = 0.98067 hPa

7. APPLICATION

Before use, the cuff pressure gauge must be tested!

Testing the cuff pressure gauge during use without connecting tube:

Seal the Luer port with a finger (see picture 1).

While the Luer port is sealed, use the manual inflation bulb to inflate the cuff pressure gauge to 100 cm H₂O (see picture 2).

The reading must remain constant for 2-3 seconds.

If the pressure drops, the device is no longer permitted to be used, and repair by the manufacturer is necessary.

When a connecting tube is used, the following leak test must be additionally performed:

Connect the connecting tube to the Luer port of the cuff pressure gauge (see picture 3).

Seal the end of the connecting tube with your finger.

While the connecting tube is sealed, use the manual inflation bulb to inflate the cuff pressure gauge to 100 cm H₂O (see picture 4).

The reading must remain constant for 2-3 seconds.

If the pressure drops, the connecting tube has a leak and must be replaced then retested.

Performing the cuff pressure check:

The cuff of the tracheostomy tube must be checked for leaks as described by the manufacturer in the instructions for use.

Before use, check the low-pressure balloon of the tracheostomy tube, and before intubation or extubation, completely evacuate the balloon by means of a syringe.

Connect the cuff pressure gauge with the connecting tube (see picture 5) or directly with the balloon supply line (see picture 6) of the tracheostomy tube.

Inflate the low-pressure balloon to a pressure of 60 to 90 cm H₂O.

Turn the red outlet valve slightly to release air until the pointer has reached the green part of the scale.

This guarantees that the balloon will press firmly against the tracheal wall.

In the case of continuous attachment to an endotracheal tube, the cuff pressure must be monitored, so that it can be adjusted to the recommended pressure range in the event of drift.

8. HYGIENE INSTRUCTIONS

Cleaning the cuff pressure gauge

Use a soft, dry cloth for cleaning.

PLEASE NOTE!

It is never permissible to immerse the device in liquids or to clean it under running water!

Disinfection can be performed by wiping, using commercial alcohol-based surface disinfectants.

The device cannot be conditioned or sterilized by machine!

Cleaning the connecting tube

The outside may be cleaned/disinfected exclusively with sterile saline solution (0.9% NaCl) followed by drying with a soft, lint-free cloth or with a suitable alcoholbased surface disinfectant compatible with the material.

PLEASE NOTE!

It is never permissible to immerse the connecting tube in liquids or to clean it under running water!

If the inside of the connecting tube is dirty or if damage is visible, the connecting tube must be properly disposed of immediately.

9. STORAGE

Store in a dry environment away from sunlight and/or heat.

10. SERVICE LIFE

The connecting tube must be replaced at the latest 7 days after the packaging has been opened or in conformity with the requirements of the healthcare standard.

11. DISPOSAL

The product must only be disposed of in accordance with the national regulations.

12. LEGAL NOTICES

The manufacturer Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for functional deficiencies, injuries, infections, and/or other complications or adverse events caused by unauthorised product alterations or improper use, care, and/or handling.

In particular, Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for damages caused by modifications to the product or by repairs, if these modifications or repairs were not carried out by the manufacturer himself.

This applies – to the extent permitted by law – both to damages to the products themselves caused thereby as well as to any consequential damages caused thereby.

When using the product beyond the period of use specified under Section 10, and/or use, treatment, maintenance (cleaning, disinfection) or storage of the product in non-compliance with the specifications laid down in these instructions for use, Andreas Fahl Medizintechnik-Vertrieb GmbH will, as far as legally permissible, be free of any liability, including liability for defects.

Should a serious incident occur in connection with this product of Andreas Fahl Medizintechnik-Vertrieb GmbH, then this must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient are domiciled.

Sale and delivery of all Andreas Fahl Medizintechnik-Vertrieb GmbH products is carried out exclusively in accordance with our General Terms and Conditions of Business (GTC) which can be obtained directly from Andreas Fahl Medizintechnik-Vertrieb GmbH.

The manufacturer reserves the right to make changes to the product at any time without notice.

MUCOPROTECT® is a trademark and brand of Andreas Fahl Medizintechnik-Vertrieb GmbH, Cologne registered in Germany and the EU member states.

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