
SPEAKING VALVE

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
PICTROGRAM LEGEND

Pictogrammes listed below you'll find on the product packaging, if applicable.

en

	With oxygen support		Keep away from sunlight
	Speaking valve (PHON)		Keep dry
	MRT suitable		Do not reuse
	Catalogue number		For reuse on a single patient
	Batch code		Non-sterile
	Use by		Medical devices
	Date of manufacture		CE marking
	Content (in pieces)		CE marking identification number of the notified body
	Consult instructions for use		Manufacturer

SPEAKING VALVE

1. FOREWORD

These instructions for use apply to **FAHL**[®] speaking valves.

These instructions for use serve to inform the physician, nursing staff and patient/user in order to ensure technically correct handling.

When used for the first time, the products must be selected, used and inserted by appropriately trained physicians or medical professionals/medical device consultants.

Please read the instructions for use carefully before using the product for the first time!

Keep them in an easily accessible place for future reference.

Please keep the package for as long as you use the product.

It contains important information on the product!

2. SAFETY INFORMATION

Tracheostomy tubes with **FAHL**[®] speaking valve must never under any circumstances be used by laryngectomised patients (patients without larynx), as this could lead to serious complications, and even to suffocation!

Under no circumstances may speaking valves be inserted directly into the tracheostoma.

Do not repair or modify the speaking valves.

Care must be taken to ensure that the valve opening is not blocked by clothing or the like in order to enable unhindered breathing (risk of respiratory distress).

The speaking valve may not be used in conjunction with a nebuliser or atomiser.

The speaking valve may not be used in conjunction with a nebuliser or atomiser.

This would invariably entail the risk of increased respiratory resistance and thus obstruction of breathing.

If the filter is clogged or if signs of clogging appear, the **FAHL**[®] speaking valves must immediately be removed and disposed of.

There is a fire hazard in connection with oxygen devices and oxygen therapies.

Do not use in the vicinity of sparks or open flames.

Do not grease connections, tubes or other accessories in order to avoid fire hazard.

Smoking during oxygen therapy is dangerous and can lead to severe fire-induced injury.

Only use lotions or ointments labelled as oxygen-compatible to avoid the risk of fire or burns.

3. INTENDED PURPOSE

The speaking valves enable spontaneously breathing tracheotomised patients with a partially or fully preserved larynx to speak without having to use their fingers if a phonation cannula in combination with a fenestrated inner cannula or a tracheostomy patch is present.

HUMIDOPHONE[®]

It filters the breathing air, thereby reducing intrusion of particles into the patient's airways.

O₂

Speaking valves with an O₂ connection also allow for the supply of additional oxygen via the integrated O₂ connection.

4. INDICATION

Speaking valves are indicated for spontaneously breathing tracheotomised patients with a partially or fully preserved larynx.

5. CONTRAINDICATION

The **FAHL**[®] speaking valve must only be used by patients with partially or fully preserved larynx!

Use in patients with strong exsiccosis (dehydration), or excessive secretion and water accumulation in the lungs and airways (left-ventricular heart failure, pulmonary oedema), or strong air leakage flow (exhaled air not flowing through the filter cassette) is not permissible.

The **FAHL**[®] speaking valve must only be used by awake and alert patients with a spontaneous breathing!

The latter does not apply if the patient is in an institution under nursing supervision for 24 hours or in home care.

The **FAHL**[®] speaking valve must not be used in case of increased secretion formation (hypersecretion).

The **FAHL**[®] speaking valve must not be used in case of pronounced vocal cord paralysis (recurrent laryngeal nerve paralysis).

6. COMPLICATIONS AND ADVERSE REACTIONS

The **FAHL**[®] must be removed before going to sleep to prevent potential obstruction of the airways!

7. PRODUCT DESCRIPTION

When breathing in, the valve membrane of the speaking valve opens.

After breathing in, the valve closes, thereby enabling the user to speak.

The pressure rise during expiration supports the valve closure process.

CAUTION!

The speaking valve is intended for use only and exclusively in combination with a phonation cannula.

COMBIPHON®

COMBIPHON® Slim: The plastic housing is equipped on the patient side with a central opening with 15 mm internal diameter, thereby allowing for connection to a 15 mm standard connector.

COMBIPHON®: The functional multi-adapter features both a central opening for attachment onto a tracheostomy tube with 15 mm connector, as well as a 22 mm combi-adapter.

COMBIPHON® O₂ & COMBIPHON® Slim O₂: The O₂ speaking valve also has an oxygen connecting piece on the side of the housing.

This allows the connection of a connecting tube and thus the combination with an oxygen source.

FAHL® BIESALSKI SPEAKING VALVE

The plastic housing is equipped on the patient side with a central opening with 15 mm internal diameter, thereby allowing for connection to a 15 mm standard connector.

HUMIDOPHONE®

The **HUMIDOPHONE®** is a speaking valve as well as a particle filter, consisting of a foam filter and a plastic housing as well as of a lid with valve function.

The plastic housing is equipped on the patient side with a standardized opening, thereby allowing for connection to a fenestrated/sieved tracheostomy tube.

HUMIDOPHONE® - On the patient side, the plastic housing is equipped with a central opening with 22 mm external diameter, thereby allowing for connection to a 22 mm receptacle system.

HUMIDOPHONE® TYP 15 - On the patient side, the plastic housing is equipped with a central opening with 15 mm internal diameter, thereby allowing for connection to a 15 mm standard connector.

HUMIDOPHONE® PLUS: The functional multi-adapter features both a central opening for attachment onto a tracheostomy tube with 15 mm connector, as well as a 22 mm combi-adapter.

HUMIDOPHONE® PLUS O₂: The functional multi-adapter features both a central opening for attachment onto a tracheostomy tube with 15 mm connector, as well as a 22 mm combi-adapter.

An oxygen connecting piece fitted to the side of the housing (5 mm plug-in connection) enables a connecting tube to be connected up, thereby allowing combination with an oxygen source.

LARYVOX® HANDS-FREE VALVE

On the patient side, the plastic housing is equipped with a central opening with 22 mm external diameter, thereby allowing for connection to a 22 mm receptacle system.

The strap running across the housing can be used for manual fine tuning of the diaphragm setting and as a spacer for textiles

The **NORMAL RESISTANCE** version has a blue valve diaphragm, which is somewhat thicker and heavier than the uncoloured valve diaphragm in the **LOW RESISTANCE** version

In the **LOW RESISTANCE** version, only low blowing pressure is required to close the speaking valve.

PHON VALVE

The plastic housing is equipped on the patient side with a central opening with 15 mm internal diameter, thereby allowing for connection to a 15 mm standard connector.

TRACHLINE® PHON VALVE: It is equipped on the patient side with a bayonet closure and thereby ensures connection to a **TRACHLINE®** tracheostomy tube.

8. APPLICATION

Check whether the packaging and the product are in a faultless condition.

If a product is damaged, it must be discarded immediately.

Check the use-by or expiry date.

Do not use the product after this date.

CAUTION!

Before using a speaking valve, the cuff of the tracheostomy tube must be completely deflated and void of air.

The speaking valve should be removed if respiratory resistance is too high.

WARNING!

Before using the speaking valve, attach the enclosed warning label to the blocking tube of the tracheostomy tube!

LARYVOX® HANDS-FREE VALVE

You should adjust the speaking valve before attaching it and set the opening thickness of the valve disc by turning it.

The further the valve is opened, the more air can be inhaled.

However, the exhalation pressure necessary to close the valve while speaking is then also greater.

For normal respiration and normal physical activity, it is therefore recommended that you set the valve such that sufficient air can flow in, but without opening the valve completely.

It is recommended that you open the valve to the maximum only for strenuous physical and sporting activities, in which your breathing becomes deeper.

After you have adjusted the valve disc to the desired opening thickness, insert the speaking valve together with the insertion ring, which can be seen on the small circumferential retaining bead, into the 22 mm adapter of the cannula.

HUMIDOPHONE®

Use with a 15 mm connector: mount the **HUMIDOPHONE®** with the central 15 mm opening by means of light pressure on the 15 mm standard connector of the tracheostomy tube (see picture 1).

Use with a 22 mm connector: mount the **HUMIDOPHONE®** by means of light pressure on the housing ring of the tracheostomy tube (see picture 2), until the housing snaps in audibly through the outer fastening rim.

Remove the **HUMIDOPHONE®** by holding the tracheostomy tube firmly with two fingers in order to extract the filter cassette by tugging lightly (see picture 3/4).

CAUTION!

To avoid injury, it is advisable to first connect the speaking valve to the tracheostomy tube and only then to insert the tracheostomy tube with attached speaking valve into the tracheostoma.

Otherwise, excessive pressure could be exerted on the mucous membranes via the cannula tube when the speaking valve is snapped into place.

PLEASE NOTE!

Do not use excessive force when connecting the speaking valve to the tracheostomy tube to avoid damage to the housing.

COMBIPHON® & FAHL® BIESALSKI SPEAKING VALVE & PHON VALVE

Carefully take the valve out of the packaging; clean and dry it prior to first use in accordance with the instructions (9. Hygiene Instructions).

The user must check the speaking valve for completeness, absence of damage and proper functioning each time before inserting the speaking valve again.

If the speaking valve is incomplete, damaged or no longer functioning properly, it must not be used anymore and must be disposed of immediately.

Use with a 15 mm connector: mount the **FAHL®** speaking valve with the central 15 mm opening by means of light pressure on the 15 mm standard connector of the tracheostomy tube.

Use with a 22 mm connector: mount the **FAHL®** speaking valve by means of light pressure on the housing ring of the tracheostomy tube, until the housing snaps in audibly through the outer fastening rim.

O₂:

For additional administration of oxygen, a connecting tube can be connected to the **O₂** connection. To this end, push the connecting tube onto the **O₂** connection.

9. HYGIENE INSTRUCTIONS

HUMIDOPHONE®

HUMIDOPHONE® are single-patient products and intended for single use only.

It must not be cleaned or disinfected, as use of cleaning agent can render the filter medium unusable. In addition, bacterial or fungicidal colonisation of the filter material can lead to an increased risk of infection.

Cleaning, disinfection or (re-)sterilisation as well as use longer than 24 hours can impair the safety and function of the product and are therefore not permissible!

COMBIPHON® & FAHL® BIESALSKI SPEAKING VALVE & LARYVOX® HANDS-FREE VALVE & PHON VALVE

Place the speaking valve to soak in clear, lukewarm water for approx. 10 minutes for cleaning. Then rinse with plenty of clear water and allow to dry in air.

PLEASE NOTE!

Do not use any sharp-edged cleaning utensils!

Heating, disinfection with boiling water or steam sterilisation are not permitted and lead to damage of the speaking valve.

10. STORAGE

Store in a dry environment away from sunlight and/or heat.

11. SERVICE LIFE

HUMIDOPHONE®

HUMIDOPHONE® are single-patient products and intended for single use only.

The maximum period of use / wearing period is 24 hours.

Speaking valves are even to be changed several times a day if required, for instance if the respiratory resistance has increased to an unacceptable level due to accumulation of secretions.

COMBIPHON® & FAHL® BIESALSKI SPEAKING VALVE & PHON VALVE

The maximum period of use must not exceed 29 days.

The service life of the product is influenced by many factors.

Thus, the composition of secretions, conscientious cleaning and storage as well as other aspects can be of decisive significance.

LARYVOX® HANDS-FREE VALVE

The maximum period of use must not exceed 365 days.

CAUTION!

A speaking valve showing damage of any kind, even if only minor, must be replaced immediately.

12. DISPOSAL

The product must only be disposed of in accordance with the national regulations.

13. LEGAL NOTICES

The manufacturer Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for damages (in particular for functional deficiencies, injuries, infections, and/or other complications or adverse events) caused by unauthorised product alterations, by repairs not performed by the manufacturer or by improper handling, care (cleaning/disinfection) and/or storage of the products in violation of the provisions of these instructions for use.

This applies – to the extent permitted by law – both to damages to the products themselves caused thereby as well as to any consequential damages caused thereby.

When using the product beyond the period of use specified under Section 10, and/or use, treatment, maintenance (cleaning, disinfection) or storage of the product in non-compliance with the specifications laid down in these instructions for use, Andreas Fahl Medizintechnik-Vertrieb GmbH will, as far as legally permissible, be free of any liability, including liability for defects.

Should a serious incident occur in connection with this product of Andreas Fahl Medizintechnik-Vertrieb GmbH, then this must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient are domiciled.

Sale and delivery of all Andreas Fahl Medizintechnik-Vertrieb GmbH products is carried out exclusively in accordance with our General Terms and Conditions of Business (GTC) which can be obtained directly from Andreas Fahl Medizintechnik-Vertrieb GmbH.


The manufacturer reserves the right to make changes to the product at any time without notice.

COMBIPHON®, **FAHL®**, **HUMIDOPHONE®**, **LARYVOX®** & **TRACHLINE®** is a trademark and brand of Andreas Fahl Medizintechnik-Vertrieb GmbH, Cologne registered in Germany and the EU member states.

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