



- ▶ Five types of fenestrated inner cannulas for DURATWIX® tracheostomy tubes, each available in two lengths:
  1. With integrated speaking valve (silicone) to avoid dislocation of speaking valve (PHON ICF)  
A cough cap releases overpressure (REF 19841 and REF 19842 short)
  2. With rotating 15 mm swivel connector (ICFV) to buffer tension and torque forces exerted by the ventilation equipment (REF 19871 and REF 19872 short)
  3. With fixed 15 mm connector (ICFU) for attachment of an HME or a speaking valve (REF 19861 and REF 19862 short)
  4. With flat profile (ICF) if no connection to any peripheral device (ventilation equipment, HME or speaking valve) is desired (REF 19851 and REF 19852 short)
  5. With 22 mm adapter (ICFK) to connect “flat profile” speaking valves or HMEs with 22 mm adapter (REF 19881 and REF 19882 short)
- ▶ Double colour coding:  
Connectors and/or connector base: turquoise colour (white colour for non-fenestrated tubes)  
Tube shaft: white colour for standard length variant; turquoise colour for short variant
- ▶ Size to be selected according to the size of the DURATWIX® outer cannula

**STERILE EO** Sterilized using ethylene oxide

Do not resterilize

Keep away from sunlight

Single patient use

Consult instructions for use

**Rx only** CAUTION: United States Federal law restricts this device to sale, distribution and use by or on order of a physician or a licensed practitioner (in the US only).

**CE 0482** CE marking identification number of the notified body

**Product made in Germany**

## Order information:

<b>DURATWIX® PHON ICF</b> , SIZE 7-10, pack of 2		<b>REF 19841</b>
<b>DURATWIX® PHON ICF short</b> , SIZE 7-10, pack of 2		<b>REF 19842</b>
<b>DURATWIX® ICFV</b> , SIZE 6-10, pack of 3		<b>REF 19871</b>
<b>DURATWIX® ICFV short</b> , SIZE 6-10, pack of 3		<b>REF 19872</b>
<b>DURATWIX® ICFU</b> , SIZE 6-10, pack of 3		<b>REF 19861</b>
<b>DURATWIX® ICFU short</b> , SIZE 6-10, pack of 3		<b>REF 19862</b>
<b>DURATWIX® ICF</b> , SIZE 7-10, pack of 3		<b>REF 19851</b>
<b>DURATWIX® ICF short</b> , SIZE 7-10, pack of 3		<b>REF 19852</b>
<b>DURATWIX® ICFK</b> , SIZE 7-10, pack of 3		<b>REF 19881</b>
<b>DURATWIX® ICFK short</b> , SIZE 7-10, pack of 3		<b>REF 19882</b>

When ordering, please state the size after the item number!