





# DURAPUR® ICF

Inner Cannulas, Fenestrated



- ▶ Five types of fenestrated inner cannulas for DURAPUR® tracheostomy tubes:
  1. With flat profile (ICF) if no connection to any peripheral device (ventilation equipment, HME or speaking valve) is desired (REF 23807)
  2. With fixed 15 mm-connector (ICFU) for attachment of an HME or a speaking valve (REF 23816)
  3. With integrated speaking valve (silver) to avoid dislocation of speaking valve (PHON ICF) (REF 23808)
  4. With integrated speaking valve (silver) to avoid dislocation of speaking valve (PHON ICF-O2) incl. O<sub>2</sub> connection (REF 23809)
  5. With 22 mm-adapter (ICFK) to connect “flat profile” speaking valves or HMEs with 22 mm-adapter (REF 23885)
  
- ▶ Size to be selected according to the size of the DURAPUR® outer cannula

 Keep away from sunlight  
 Single patient use  
 Consult instructions for use  
 **CE 0482** CE marking identification number of the notified body

## Order information:

<b>DURAPUR® ICF</b>	SIZE 7-12	  	<b>REF 23807</b>
<b>DURAPUR® ICFU</b>	SIZE 7-12	  	<b>REF 23816</b>
<b>DURAPUR® PHON ICF</b>	SIZE 7-12	 	<b>REF 23808</b>
<b>DURAPUR® PHON ICF-O2</b>	SIZE 7-12	  	<b>REF 23809</b>
<b>DURAPUR® ICFK</b>	SIZE 7-12	  	<b>REF 23885</b>

When ordering, please state the size after the item number!