

COMBIPHON® / COMBIPHON® O₂

REF 27131



CE

REF 27132



CE 0482

- ▶ Speaking valve - transparent, lightweight
- ▶ With silicone membrane
- ▶ Membrane closed in basic position ("closed-position valve") – opens during inspiration. Advantage: No enhanced exhalation pressure required to close the valve and to facilitate speaking. Less fatigue for the patient, better speaking quality
- ▶ With multiadapter: fits to both 15 mm connector (**UNI**) and 22 mm adapter (**KOMBI**) tracheostomy tubes or base plates
- ▶ To be used by tracheotomized patients, not to be used by laryngectomized patients
- ▶ To be used only and exclusively in combination with a phonation cannula.

SCOPE OF DELIVERY

COMBIPHON® REF 27131

1 x speaking valve
1 x Instruction for use

or

COMBIPHON® O₂ REF 27132

1 x speaking valve
1 x Instruction for use

Product	Unit	REF
COMBIPHON®	1 piece	27131
COMBIPHON® O ₂	1 piece	27132



15 mm

UNI



22 mm

KOMBI

O₂O₂

The COMBIPHON® speaking valve must only be used by patients with partially or fully preserved larynx!