



TRACHEOTEC® VARIO / TRACHEOTEC® VARIO CUFF

- ▶ Tracheostomy Tube with (REF 18920) and without CUFF (REF 18560); also available with suction line (REF 18930; on separate data sheet)
- ▶ Soft, flexible neck-flange with key tube identification data clearly marked
- ▶ With rotating 15 mm connector (VARIO) to buffer torque and tension caused by tractive forces from ventilation equipment
- ▶ With X-ray contrast stripe
- ▶ Set includes tube holder, obturator and product ID card



Order information:

TRACHEOTEC® VARIO, SIZE 3-9

TRACHEOTEC® VARIO CUFF, SIZE 3-10

When ordering, please state the size after the item number!


REF 18560

REF 18920

Size	O.D. neck flange (mm)	O.D. Tip (mm)	I.D. (mm)	Length midline (mm)	Length (mm)	REF	REF mit CUFF
3	4.7	4.7	3.0	43.7	47.2	18560-03	18920-03
4	6.0	6.0	4.0	52.7	57.1	18560-04	18920-04
5	7.3	7.3	5.0	55.7	61.1	18560-05	18920-05
6	8.7	8.7	6.0	60.7	67.1	18560-06	18920-06
7	10.0	10.0	7.0	69.6	77.0	18560-07	18920-07
8	11.3	11.3	8.0	79.5	87.9	18560-08	18920-08
9	12.7	12.7	9.0	85.5	94.9	18560-09	18920-09
10	14.0	14.0	10.0	95.0	105.3	-	18920-10

O-tube = outer cannula · O.D. neck flange = outer diameter behind the neck flange · O.D. Tip = outer diameter at the tube tip · I.D. = inner diameter at cannula tip · I-tube = inner cannula · Length = length over the outer curve · Length at mid radius = length of centre line from neck flange to tip