



# TRACHLINE® UNI CUFF

1 IC + 2 ICU

- ▶ Tracheostomy tube with CUFF for long-term use, max. service life up to 29 days
- ▶ Made of thermosensitive medical synthetic materials
- ▶ Soft, pliable, lightweight
- ▶ Transparent, inconspicuous design
- ▶ Rounded tube tip to prevent irritation of the mucous membranes
- ▶ Bayonet locks enable simple, secure attachment of the inner cannula and accessories
- ▶ Set consists of an inner cannula (IC), two inner cannulas with 15 mm connector (ICU), a 15 mm connector (UNI), a cough cap, a locking ring, an adjustable tube holder, an obturator as well as a cannula passport



## Order information:

**TRACHLINE® UNI CUFF 1 IC + 2 ICU**      Size 5-10



**REF 13800**

When ordering, please state the size after the item number!

Size	O.D. neck flange mm	I.D. O.C. mm	I.D. I.C. mm	Angel ca.	Cuff (Ø) mm	Length mm
5	10,6	6.5	5.0	115°	18	62
6	11,7	7.5	6.0	110°	22	70
7	12,8	8.5	7.0	105°	22	72
8	13,8	9.5	8.0	100°	26	80
9	14,8	10.5	9.0	95°	30	82
10	15,8	11.5	10.0	90°	32	84

O.D. neck flange = Outer diameter behind the neck flange · I.D. = Inner diameter at the tip of the cannula · O.C. = Outer cannula · I.C. = Inner cannula · Length = Length over middle curve