

Blom-Singer® SpeakFree™ HME Hands Free Valve

REF BE 1090EZ



REF BE 1090EF



- ▶ Blom-Singer® SpeakFree™ HME Hands Free Valve is a heat and moisture exchange (HME) cartridge with a single-use fully integrated adjustable hands free valve.
- ▶ Easy change between hands free or manual occlusion
- ▶ Supports a good breathability and moisture exchange
- ▶ Available in two versions:
 - **ClassicFlow (BE 1090EZ)** with normal foam
 - **EasyFlow (BE 1090EF)** with open-cell foam
- ▶ Is working as a hygroscopic filter
- ▶ Attachable crossbar prevents accidental airway obstruction.
- ▶ With 22 mm adapter (KOMBI) to connect to any conventional base plates (e.g. LARYVOX® TAPE) or silicone tubes for laryngectomees (e.g. LARYNGOTEC® PRO)
- ▶ Made with latex-free materials
- ▶ Discreet thanks to a non-reflective surface texture

SCOPE OF DELIVERY

- ▶ Instruction for use

**Blom-Singer®
SpeakFree™ HME Hands Free Valve
Classic Flow** **REF BE 1090EZ**

30 x Hands Free speech valve
5 x crossbars

or

**Blom-Singer®
SpeakFree™ HME Hands Free Valve
Easy Flow** **REF BE 1090EF**

30 x Hands Free speech valve
5 x crossbars

Product	Unit	REF
Blom-Singer® SpeakFree HME™, Classic Flow	30 pieces	BE 1090EZ
Blom-Singer® SpeakFree HME™, Easy Flow	30 pieces	BE 1090EF

Andreas Fahl Medizintechnik-Vertrieb GmbH is the exclusive distributor for Blom-Singer® products in selected European countries.
Blom-Singer® is a registered trademark of Freudenberg Medical, LLC.