


OPTIFLAUSCH®

REF 32550



REF 32551



- ▶ Tracheostomy tube holder
 - "K" variant: with adherent fastener
 - "H" variant: with hook fastener
- ▶ Assures the secure fixation of tracheostomy tubes
- ▶ Foam core surrounded by soft non-woven fabric
- ▶ Broad and soft for high wearing comfort
- ▶ Approx. 4 mm thick
- ▶ The tube holder is wrapped around the neck and the adherent fastener enables it to be easily and quickly adjusted to the desired neck size (2 parts, adjustable).
- ▶ Colour: white
- ▶ Individually packed

SCOPE OF DELIVERY

OPTIFLAUSCH® K **REF 32550**

1 x Tracheostomy Tube Holder
1 x Instruction for use

or

OPTIFLAUSCH® H **REF 32551**

1 x Tracheostomy Tube Holder
1 x Instruction for use

| Product | Length | Width | REF |
|-----------------------|------------|--------|--------------|
| OPTIFLAUSCH® K | 24 - 45 cm | 2,5 cm | 32550 |
| OPTIFLAUSCH® H | 24 - 47 cm | 2,5 cm | 32551 |



OPTIFLAUSCH®

REF 32553



REF 32555



- ▶ Tracheostomy tube holder
 - “K” variant: with adherent fastener
 - “K PED” variant: with adherent fastener, for pediatric use
- ▶ Assures the secure fixation of tracheostomy tubes
- ▶ Foam core surrounded by soft non-woven fabric
- ▶ Broad and soft for high wearing comfort
- ▶ Approx. 4 mm thick
- ▶ The tube holder is wrapped around the neck and the adherent fastener enables it to be easily and quickly adjusted to the desired neck size (2 parts, adjustable).
- ▶ With white adherent fastener and elasticated insert
- ▶ Colour: beige
- ▶ Individually packed

SCOPE OF DELIVERY

OPTIFLAUSCH® K beige REF 32553

1 x Tracheostomy Tube Holder
1 x Instruction for use

or

OPTIFLAUSCH® K PED REF 32555

1 x Tracheostomy Tube Holder
1 x Instruction for use

| Product | Length | Width | REF |
|----------------------|------------|--------|-------|
| OPTIFLAUSCH® K beige | 24 - 45 cm | 2,5 cm | 32553 |
| OPTIFLAUSCH® K PED | 16 - 23 cm | 2,0 cm | 32555 |

