

TRACHFLOW® SHUNT

sterile

REF 68800



REF 60700



- ▶ Allows a targeted aspiration of shunt valves / voice prosthesis
- ▶ Made of soft, flexible and transparent plastic
- ▶ Colour-coded connectors
- ▶ Length 10 cm
- ▶ To be used in conjunction with the TRACHFLOW® Fingertip to control suction strength.
- ▶ The TRACHFLOW® Fingertip is used to connect the suction tube of the suction device to the sterile stoma catheter.
- ▶ The conical design allows one size to fit all sizes of stoma catheters

SCOPE OF DELIVERY


- ▶ Manual

TRACHFLOW® SHUNT **REF 68800**
1 x Suction catheter, sterile

When ordering, please state the size after the item number!

Optional accessories:
TRACHFLOW® FINGERTIP **REF 60700**

INFORMATION



Article	Size	Length	Color	Unit	REF
TRACHFLOW® SHUNT	CH 08	10 cm	Blue	30	68800-08
TRACHFLOW® SHUNT	CH 10	10 cm	Black	30	68800-10

CE 0482 STERILE EO