

2-KAM® SLIT STERILE

sterile

REF 30606



- ▶ Sterile double-chamber tracheal compress
- ▶ The 2-KAM® compress consists of two compress layers of different thicknesses, which are fixed on each other and made of different non-woven fabrics
- ▶ The top of compress body is reinforced and can be wiped without linting, and absorbs leaking secretion quickly
- ▶ The back of the compress is particularly soft and keeps the skin dry and warm
- ▶ Tracheal compresses are used to absorb tracheal secretions as well as to cushion the tracheal neck flange against the skin
- ▶ It is sterile and can therefore be used after an operation
- ▶ The 2-KAM® **SLIT** variant serves for a quick change of compresses while the tracheostomy tube is in situ
- ▶ For single use only

SCOPE OF DELIVERY

- ▶ Manual

2-KAM® SLIT STERILE**REF 30606**

100 x tracheal compresses, sterile

Product	Size	Unit	REF
2-KAM® SLIT STERILE	9,0 x 9,8 cm	100	30606

CE 0482 STERILE EO