

OPTIFLUID® STOMA OIL

- ▶ Colourless, neutral oil to lubricate tracheostomy tubes
- ▶ Non-irritant to the mucosa
- ▶ Facilitates the insertion of tracheostomy tubes into the tracheostoma
- ▶ Also to be used for tracheostomy tubes with inner cannulas: Inner cannulas may be lubricated to facilitate their insertion into the outer cannula
- ▶ Economic use: only a few drops needed for proper lubrication of inner and/or outer cannula
- ▶ PET-bottle with dropper to allow sparing and exact dosage of the oil
- ▶ Also available as stoma oil towel:
no dripping of the oil, assurance of the application of the proper quantity,
individually packed in one-each sachets; ideal for travelling



Order Information:

OPTIFLUID® Stoma oil, Bottle 25 ml

REF 31525

OPTIFLUID® Stoma oil towel, Pack of 30 Stück

REF 31550