

# BLOM-SINGER® ELECTROLARYNX

Digital Speech Aid

- ▶ Artificial larynx which is used to restore speech in post-laryngectomy patients
- ▶ Allows voice restoration after surgery
- ▶ Easy to learn and operate: The pitch and volume of the voice are simply adjustable via two dial switches to meet the individual needs of the patient
- ▶ From 40 Hz to 170 Hz
- ▶ Operated by standard 9V block batteries; does not require a charging device so it is not dependent on a power supply.
- ▶ Complemented by an oral adapter to introduce sound directly to the back of the oral cavity when the sound does not transmit effectively through the tissues of the side of the throat following surgery and/or radiation



Consult instructions for use

## Order information:

**Blom-Singer® ElectroLarynx Digital Speech Aid**

**REF EL1000**