


LARYVOX® EXTRA HME HIGHFLOW STARTER SET

Set consisting of:

- ▶ 3x LARYVOX® TAPE STANDARD OVAL
base plate with water-repellent effect, with 22 mm adapter (KOMBI)
- ▶ 3X LARYVOX® TAPE FLEXIBLE OVAL
base plate of soft, textile-type material with high adhesive force, with 22 mm adapter (KOMBI)
- ▶ 3X LARYVOX® TAPE HYDROSOFT OVAL
base plate of flexible and skin-compatible material based on hydrocolloids, with 22 mm adapter (KOMBI)
- ▶ 3X LARYVOX® TAPE COMFORT OVAL
base plate of thin, breathable skin-coloured material, with 22 mm adapter (KOMBI)
- ▶ 3x LARYVOX® TAPE EXTRA FINE rectangular
base plate particularly suitable for sensitive skin, with 22 mm adapter (KOMBI)
- ▶ 15x OPTICLEAR®
Special wipe for removal of all types of skin adhesive, particularly for use around the tracheostoma;
Soft cleaning fleece with an alcohol-based cleaning lotion, adhesive residues are dissolved and removed from the skin gently and atraumatically
- ▶ 15x OPTIGARD®
Special wipe for preparatory skin care in the use of adhesive plasters/base plates and skin adhesives
- ▶ 15x LARYVOX® EXTRA HME HighFlow
Beige HME with 22 mm adapter (KOMBI) for connection to all standard base plates;
HighFlow model: reduced respiratory resistance under high physical stress
- ▶ 1x LARYVOX® SECUTRACH® shower guard
Shower guard; consists of a plastic housing and serves to protect the tracheostoma while showering and washing hair.
The shower guard prevents water from entering the tracheostoma and is designed for use with a 22 mm connector ring.



Please refer to the pictograms in the respective instructions of each product.

Order information:

LARYVOX® EXTRA HME HIGHFLOW STARTER SET

REF 48840