


# NECKFIX®

REF 32620



REF 32610



REF 32612



REF 32622



- ▶ Tracheostomy tube holder
  - “K” variant: with adherent fastener - white
  - “Blue K” variant: with adherent fastener - blue
  - “Blue PED K” variant: with adherent fastener, for pediatric use – blue
  - “Beige K PED” variante: with adherent fastener, for pediatric use – beige
  - “Beige K” variante: with adherent fastener – beige
  - “Beige K XL” variante: with adherent fastener – beige
- ▶ Assures the secure fixation of tracheostomy tubes
- ▶ Foam core surrounded by soft non-woven fabric
- ▶ Economic tube holder for routine use
- ▶ Slim design, one part
- ▶ Individually packed

## SCOPE OF DELIVERY

▶ Manual	
<b>NECKFIX® K</b>	<b>REF 32620</b>
1 x Tracheostomy Tube Holder	
oder	
<b>NECKFIX® Blue K</b>	<b>REF 3261x</b>
1 x Tracheostomy Tube Holder	
oder	
<b>NECKFIX® Beige K</b>	<b>REF 3262x</b>
1 x Tracheostomy Tube Holder	

Product	Lenght	Width	REF
NECKFIX® K	45 cm	2,5 cm	32620
NECKFIX® Blue K	45 cm	3,5 cm	32610
NECKFIX® Blue PED K	34 cm		32612
NECKFIX® Beige K PED	34 cm		32621
NECKFIX® Beige K	46 cm		32622
NECKFIX® Beige K XL	50 cm		32623

