

LARYVOX® TAPE EXTRA FINE

Base Plate

- ▶ Very good adhesive properties
- ▶ Non-invasive
- ▶ Transparent and extremely thin base plate
- ▶ Water-repellent material
- ▶ Covers a large skin area around the tracheostoma
- ▶ Adapts optimally to the movements of the neck region
- ▶ Discreet appearance
- ▶ No friction inside the tracheostoma
- ▶ With 22 mm housing to be combined with any HME and speaking valve with 22 mm adapter (KOMBI)
- ▶ Can be combined with special silicone tubes (LARYNGOTEC® KOMBI CLIP or LARYNGOTEC® KOMBI CLIP LINGO) if additional tracheal stabilisation is needed
- ▶ Shape: rectangular



Consult instructions for use



CE marking



CE-Kennzeichnung

Product made in Germany

Order information:

LARYVOX® Tape EXTRA FINE, pack of 15

REF 48400