

HUMIDOTRACH® / HUMIDOTRACH® ISO



- ▶ HME filter with oxygen supply port
- ▶ 15 mm standard connector fits to all tracheostomy tubes with 15 mm connector (UNI)
- ▶ Open pore polyurethane foam filter
- ▶ Grid shaped housing to facilitate a maximum amount of air exchange
- ▶ With overpressure valve to allow unimpaired coughing
- ▶ Individually packed in convenient peel-packs
- ▶ **HUMIDOTRACH®**: An oxygen connecting piece fitted to the side of the housing (5 mm plug-in connection) enables a connecting tube to be connected up, thereby allowing combination with an oxygen source
- ▶ Oxygen connecting tube separat available (REF 46320)
- ▶ **HUMIDOTRACH® ISO**: An oxygen connecting piece fitted to the side of the housing (6 mm ISO plug-in connection / nipple) enables a connecting tube to be connected up, thereby allowing combination with an oxygen source
- ▶ Standardised for all ISO oxygen connection tubes

Order information:

HUMIDOTRACH®	Pack of 30	REF 46450
HUMIDOTRACH® ISO	Pack of 30	REF 46455
Oxygen connecting tube	Pack of 50	REF 46320