

HUMIDOSTOM®/ HUMIDOSTOM® O₂/ HUMIDOSTOM® O₂ ISO

REF 46160



REF 46150



REF 46155



- ▶ HME filter (Heat-Moisture Exchanger)
- ▶ 15 mm standard connector (UNI) fits to all tracheostomy tubes with 15 mm connector
- ▶ Open pore polyurethane foam filter
- ▶ **HUMIDOSTOM® O₂ (REF 46150)**: An oxygen connecting piece fitted to the side of the housing (5 mm plug-in connection) enables a connecting tube to be connected up, thereby allowing combination with an oxygen source
- ▶ **HUMIDOSTOM® O₂ ISO (REF 46155)**: An oxygen connecting piece fitted to the side of the housing (6 mm ISO plug-in connection / nipple) enables a connecting tube to be connected up, thereby allowing combination with an oxygen source
- ▶ Lightweight and inconspicuous design

SCOPE OF DELIVERY

- ▶ Instruction for use

HUMIDOSTOM® REF 46160

1 x HME filter cassettes

or

HUMIDOSTOM® O₂ REF 46150

1 x HME filter cassettes

or


HUMIDOSTOM® O₂ ISO REF 46155

1 x HME filter cassettes



15 mm

UNI

O₂O₂

Product	Unit	REF
HUMIDOSTOM®	30 pieces, individually packed	46160
HUMIDOSTOM® O₂	30 pieces, individually packed	46150
HUMIDOSTOM® O₂ ISO	30 pieces, individually packed	46155

CE 0482