

TRACHEOSTOMA CATHETER METRIC S + L

sterile

REF 68000/68900



REF 68100



REF 60700



- ▶ Sterile tracheostoma catheters for suctioning tracheal secretions
- ▶ Made of soft, flexible and transparent plastic
- ▶ With central opening and four lateral openings ("eyes") at the tip to avoid adherence to the trachea
- ▶ With grading in 1 cm increments to facilitate control of the insertion depth
- ▶ Wide range of sizes (CH 6-18)
- ▶ Colour coded connectors
- ▶ Available in two different lengths: 25 cm and 50 cm (the short variant is particularly suited for tracheal suctioning)
- ▶ To be used with the TRACHFLOW® Fingertip, to control the suction power. The Fingertip is connected between the suction device and the sterile stoma catheter. One size of the TRACHFLOW® Fingertip fits all sizes of stoma catheters

SCOPE OF DELIVERY

▶ Manual

TRACHEOSTOMA CATHETER METRIC S REF 68100
1 x Suction catheter, sterile, atraumatic

or

STOMA CATHETER METRIC L REF 68000
1 x Suction catheter, sterile, atraumatic

or

ABSAUGKATHETER STANDARD METRIC L REF 68900
1 x Suction catheter, sterile, atraumatic

Optional accessories:

TRACHFLOW® FINGERTIP REF 60700

INFORMATION

06 08 10 12 14 16 18

CE 0482 STERILE EO

Article	Size	Length	Shore	Unit	REF
TRACHEOSTOMA CATHETER METRIC S	CH 08-16	25 cm	A64	30	68100-xx
STOMA CATHETER METRIC L	CH 06-18	50 cm	A64	30	68000-xx
ABSAUGKATHETER STANDARD METRIC L	CH 06-18	50 cm	A78	100	68900-xx