

BLOM-SINGER® DUCKBILL VOICE PROSTHESIS

- ▶ The Blom-Singer® Duckbill voice prosthesis can be changed by the patients themselves after receiving professional instruction and briefing
- ▶ For daily cleaning, the voice prosthesis is removed from the shunt channel / fistula channel and then replaced
- ▶ Fixation of the security thread outside of the tracheostoma is absolutely necessary
- ▶ The voice prosthesis, which has a uniform shaft thickness of 16 Fr., is available in sizes 6 - 25



Single patient use



Consult instructions for use

Order information:

Blom-Singer® Duckbill - 16 Fr.

SIZE 6	REF DB16-006
SIZE 8	REF DB16-008
SIZE 10	REF BE6000
SIZE 12	REF DB16-012
SIZE 14	REF BE6001
SIZE 18	REF BE6002
SIZE 25	REF BE6004