

# LARYVOX® UNIEK ULTRA LOW RESISTANCE

sterile

REF 25010



REF 25020



REF 25925



REF 25800 (optional)

Please scan  
for the videos

CE 0482

- ▶ A choice of 8 length sizes (5 - 11 mm / 13 mm) and 2 diameters (21 Fr. and 24 Fr.) is available
- ▶ Made of medical silicone and radiocontrast barium sulphate to allow for clear visualisation in case of aspiration / ingestion
- ▶ The LARYVOX® UNIEK U.L.R. voice prosthesis is suitable for MRI examinations
- ▶ Placement is possible anterograde with the LARYVOX® INSERTER UNIEK (REF 25925)
- ▶ The LARYVOX® UNIEK U.L.R. is supplied packaged sterile
- ▶ The variant of the 24 Fr. LARYVOX® UNIEK U.L.R. voice prosthesis is suitable for specific use in patients with a distended shunt canal
- ▶ The LARYVOX® PLUG (REF 25800) offers a temporary solution in case of possible transprosthetic leakages
- ▶ Video tutorials can be viewed on YouTube for a better understanding, for example, for cleaning and changing

## SCOPE OF DELIVERY

- ▶ Manual, voice prosthesis pass

### LARYVOX® UNIEK ULTRA LOW RESISTANCE 21 Fr. REF 25010

1 x voice prosthesis  
1 x LARYVOX® INSERTER UNIEK  
1 x LARYVOX® DUO BRUSH

or

### LARYVOX® UNIEK ULTRA LOW RESISTANCE 24 Fr. REF 25020

1 x voice prosthesis  
1 x LARYVOX® INSERTER UNIEK  
1 x LARYVOX® DUO BRUSH

Optional available:

### LARYVOX® PLUG REF 25800

The correct selection and change of the voice prosthesis is done by the physician.

Size	Shaft thickness Ø mm 21 FR.	Shaft thickness Ø mm 24 FR.	Length mm	REF 21 FR.	REF 24 FR.
05	7,0	8,0	5,0	25010-05	25020-05
06	7,0	8,0	6,0	25010-06	25020-06
07	7,0	8,0	7,0	25010-07	25020-07
08	7,0	8,0	8,0	25010-08	25020-08
09	7,0	8,0	9,0	25010-09	25020-09
10	7,0	8,0	10,0	25010-10	25020-10
11	7,0	8,0	11,0	25010-11	25020-11
13	7,0	8,0	13,0	25010-13	25020-13