


TRACHEOTEC® SILC VARIO CUFF

- ▶ Tracheostomy tube with High-Pressure Cuff (CUFF)
- ▶ Soft, flexible neck-flange with key tube identification data clearly marked
- ▶ With rotating 15 mm connector (VARIO) to buffer torque and tension caused by tractive forces from ventilation equipment
- ▶ The high pressure cuff (CUFF) is not preformed and is close against the tube in the relaxed state, allowing easy and atraumatic insertion of the cannula
- ▶ Rounded cannula tip to prevent irritation of the mucosa in the trachea
- ▶ The valve of the control balloon is colour-coded and varies depending on the cannula size (see table below)
- ▶ Set includes tube holder, obturator and product ID card



Les informations de commande :

TRACHEOTEC® SILC VARIO CUFF, SIZE 6-9

When ordering, please state the size after the item number!



REF 18850

Size	O.D. Neck flange mm	O.D. Tip mm	Length mm	Colour code	REF
5.0	7.3	7.3	59	violet	18850-05
6.0	8.7	8.7	70	red	18850-06
7.0	10.0	10.0	75	green	18850-07
8.0	11.0	11.0	81	blue	18850-08
9.0	12.3	12.3	92	white	18850-09

O.D. neck flange = outer diameter behind the neck flange · O.D. Tip = outer diameter at the tube tip · O-tube = outer cannula · Length = length over the outer curve