


FAHL® DECANNULATION TAPE FLEXIBLE

REF 48001-01



REF 48001-02



REF 48001-03



- ▶ Decannulation tape for the weaning from a tracheostomy tube for tracheostomised patients
- ▶ Perfect fit to cover the tracheostoma
- ▶ The large contact surface and the cuff seal the tracheostoma securely and prevent air from escaping during speaking
- ▶ Breathable and thin material for maximum wearing comfort and best adhesive properties
- ▶ Available in three different variants in white colour:
 - LARGE: 200 mm x 155 mm
 - SMALL: 165 mm x 135 mm
 - OVAL: 140 mm x 100 mm
- ▶ The FAHL® DECANNULATION TAPE may only be used for tracheostomised patients and is contraindicated for laryngectomised patients!

SCOPE OF DELIVERY

- ▶ Syringe (vol. 20 ml)
- ▶ Instruction for use

FAHL® DECANNULATION TAPE FLEXIBLE REF 48001

5 x Decannulation Tapes with Cuff, white

When ordering, please state the size after the item number!

Product	Unit	REF
FAHL® DECANNULATION TAPE FLEXIBLE, LARGE	5	48001-01
FAHL® DECANNULATION TAPE FLEXIBLE, SMALL	5	48001-02
FAHL® DECANNULATION TAPE FLEXIBLE, OVAL	5	48001-03

