








LARYVOX® PLUG

Sealing plug



- ▶ The LARYVOX® PLUG sealing plug serves for the short-term use with a defective shunt valve to prevent transprosthetic fluid leakage (centrally through the shunt valve)
- ▶ The sealing plug should only be used temporarily until the shunt valve has been professionally inspected and replaced if necessary
- ▶ Use of the LARYVOX® PLUG prevents leakage through the voice prosthesis while eating and drinking
- ▶ Talking is not possible while the plug is inserted
- ▶ The safety strap prevents aspiration of the sealing plug
- ▶ The LARYVOX® PLUG is supplied non-sterile and with an insertion rod
- ▶ Available in two sizes for different shunt valves



Order information:

LARYVOX® PLUG

for all common shunt valves (20 - 22.5 Fr)

REF 25800

LARYVOX® PLUG 4 mm

for all common shunt valves (4 mm)

REF 25810