


TRACHEOTEX® SCARF PED

for children

- ▶ Covers and protects the tracheostoma for laryngectomised and tracheotomised patients
- ▶ It is tied around the neck with a strap (with velcro closure)
- ▶ Notes on the material quality: normal to particularly air permeability, made of cotton jersey; soft quality, very comfortable on the skin
- ▶ The centrally incorporated ventilation insert, made of multiple layers of fine mesh fabric, warms, humidifies and filters the inspired air
- ▶ The neck width can be adjusted individually to any size between 22 cm and 29 cm
- ▶ Washable (hand wash) and reusable
- ▶ Different colours (colour sample card available as a separate brochure)



Order information:

TRACHEOTEX® Scarf PED, 22 - 29 cm

Please specify color number on the color sample card when ordering!

REF 42010