



- ▶ To prevent ingress of water into the trachea during showering
- ▶ Air openings with large diameters located in the lower area of the cap to prevent ingress of water while providing sufficient air to allow breathing without anxiety
- ▶ Sides of the cap have a thin protective edge which fits seamlessly to the skin, thereby optimizing the protection against the ingress of water
- ▶ Flexible and soft material allows the SECUTRACH® shower guard to adapt comfortably even to difficult and uneven neck anatomies
- ▶ Integrated struts enable the SECUTRACH® shower guard to retain its shape and to prevent it from collapsing
- ▶ Adjustable to different neck sizes



Single patient use only



Consult instructions for use



CE marking

Order information:

SECUTRACH® shower guard

REF 47000