

TRACHFLOW® CONT

REF 65940




REF 65920



- ▶ The TRACHFLOW® CONT is a practical container for storing suction catheters in a collected and protected manner.
- ▶ The spacious opening of the lid and the transparent body allow quick removal of the desired catheter.
- ▶ The lid can be unscrewed for quick and easy cleaning.
- ▶ Different clips on the lid, for example, allow fixation to a standard rail or individual fastening by tape. Tubing systems may also be applied in an orderly manner.
- ▶ The resilient materials ensure a long-lasting service life of the holster.
- ▶ The TRACHFLOW® CONT is available in two different sizes.

SCOPE OF DELIVERY

TRACHFLOW® CONT, short	REF 65920
1 x Catheter holster, 30 cm	
1 x Screw lid	
or	
TRACHFLOW® CONT, large	REF 65940
1 x Catheter holster, 50 cm	
1 x Screw lid	



Product	Size	Unit	REF
TRACHFLOW® CONT, short	30 cm	individually packed	65920
TRACHFLOW® CONT, large	50 cm	individually packed	65940