


SUPRAPHIX®

REF 32300



REF 32500



- ▶ Tracheostomy tube holder – 4 variants:
 - “H” variant: with hook fastener
 - “H PED” variant: for pediatric use, with hook fastener
 - “K” variant: with adherent fastener
 - “K PED” variant: for pediatric use, with adherent fastener
- ▶ Assures the secure fixation of tracheostomy tubes
- ▶ Foam core surrounded by soft non-woven fabric
- ▶ Broad and soft for high wearing comfort
- ▶ Approx. 4 mm thick and 2.5 cm wide.
- ▶ The tube holder is wrapped around the neck and the adherent fastener enables it to be easily and quickly adjusted to the desired neck size (2 parts, adjustable).
- ▶ Colour: white
- ▶ Individually packed

SCOPE OF DELIVERY

- ▶ Manual

SUPRAPHIX® xx H **REF 323xx**
1 x Tracheostomy Tube Holder

or

SUPRAPHIX® xx K **REF 325xx**
1 x Tracheostomy Tube Holder

Product	Length	REF
SUPRAPHIX® 11 H PED	15 - 26 cm	32308
SUPRAPHIX® 13 H PED	17 - 28 cm	32310
SUPRAPHIX® 15 H PED	19 - 30 cm	32311
SUPRAPHIX® 18 H	22 - 41 cm	32301
SUPRAPHIX® 25 H	35 - 42 cm	32300
SUPRAPHIX® 30 H	34 - 54 cm	32305
SUPRAPHIX® 36 H	40 - 58 cm	32303
SUPRAPHIX® 11 K PED	18 - 25 cm	32509
SUPRAPHIX® 13 K PED	20 - 27 cm	32510
SUPRAPHIX® 15 K PED	22 - 29 cm	32512
SUPRAPHIX® 18 K	20 - 40 cm	32502
SUPRAPHIX® 25 K	32 - 48 cm	32500
SUPRAPHIX® 30 K	37 - 53 cm	32501
SUPRAPHIX® 36 K	43 - 57 cm	32503

