

FAHL® DECANNULATION-SET

Contents:

▶ LARYVOX® TAPE EXTRA FINE

15 pieces, transparent and extremely thin base plate with excellent adhesive properties for optimal adaption to movements of the neck

▶ FAHL® MULTI PLUG

1 piece, the FAHL® MULTI PLUG is compatible with all common 22 mm, e.g. LARYVOX® Tape base plate and 15 mm-mounting systems, e.g. SPIRAFLEX® UNI

▶ COMBIPHON® O₂

1 piece, speaking valve - closed position valve with silicone membrane - with multiadapter to allow connection to 15 mm connectors (UNI) and 22 mm adapters (KOMBI).



For all components of this kit separate and detailed data sheets are available!

Order information:

FAHL® Decannulation-Set

REF 48005