

BLOM-SINGER® ADVANTAGE INDWELLING VOICE PROSTHESIS

- ▶ Blom-Singer® ADVANTAGE Indwelling voice prosthesis with low-pressure valve
- ▶ In this version, silver oxide-containing silicone (7%) was used in the valve section in order to reduce candida-related early valve replacements
- ▶ The Blom-Singer® ADVANTAGE Indwelling is available in two different variants: hard- and soft valve
- ▶ The Blom-Singer® ADVANTAGE Indwelling - Hard Valve Assembly voice prosthesis is additionally furnished with a titanium ring for reinforcement of the shaft and is available as 20 Fr variant in size 4 - 14
- ▶ The Blom-Singer® ADVANTAGE Indwelling - Soft Valve variant has no titanium ring and is available as 16 Fr variant in size 4 - 14
- ▶ With the help of the LARYVOX® INSERTER the voice prosthesis can be placed anterograde quickly and easily:

LARYVOX® INSERTER (REF 25900) - for different types of voice prosthesis

LARYVOX® INSERTER EXPERT (REF 25910) - for 16 Fr. voice prosthesis (or DUAL VALVE™ 20 Fr.)

LARYVOX® INSERTER EXPERT SMALL (REF 25915) - for 4 mm voice prosthesis

Application with gel capsule technique:

LARYVOX® INSERTER SPECIAL - (REF 25920) - for difficult tracheostoma



Hard Valve



Soft Valve

Order information:

Blom-Singer® ADVANTAGE Indwelling Hard Valve Assembly Voice Prosthesis	20 Fr.	SIZE 4 - 14	REF AD20XX-NS
---	---------------	--------------------	----------------------

Blom-Singer® ADVANTAGE Indwelling Soft Valve Assembly Voice Prosthesis	16 Fr.	SIZE 4 - 14	REF AS16XX-NS
---	---------------	--------------------	----------------------

LARYVOX® INSERTER	REF 25900
--------------------------	------------------

LARYVOX® INSERTER EXPERT	REF 25910
---------------------------------	------------------

LARYVOX® INSERTER EXPERT SMALL	REF 25915
---------------------------------------	------------------

LARYVOX® INSERTER SPECIAL	REF 25920
----------------------------------	------------------

	Single patient use
	Consult instructions for use