

# LARYNGOTEC® KOMBI / LARYNGOTEC® KOMBI LINGO

- ▶ Tracheostomy tube for long-term use, can be used up to 3-6 months
- ▶ Medical grade silicone: soft, pliable, lightweight
- ▶ Neck flange perfectly adapted to neck anatomy, with lateral openings for fixation of tube holder
- ▶ Transparent, inconspicuous design
- ▶ With 22 mm adapter (KOMBI)
- ▶ Each size in three different lengths
- ▶ With (LARYNGOTEC® KOMBI LINGO REF 14720) or without (LARYNGOTEC® KOMBI REF 14920) sieve fenestration
- ▶ Custom-made sieve-fenestration possible
- ▶ Including broad adjustable OPTIFLAUSCH® tube holder and a product ID card



## Order information:

LARYNGOTEC® KOMBI, SIZE 7-13



REF 14920

LARYNGOTEC® KOMBI LINGO, SIZE 7-13



REF 14720

When ordering, please state the size after the item number!

Size	O.D. tip mm	O.D. neck flange mm	I.D. mm	Length mm	Length middle radius mm	REF	REF LINGO
7	10.2	10.2	7.2	44	38	14920-0744	14720-0744
7	10.2	10.2	7.2	62	54	14920-0762	14720-0762
7	10.2	10.2	7.2	65	56	14920-0765	14720-0765
8	11.2	11.2	8.2	44	38	14920-0844	14720-0844
8	11.2	11.2	8.2	62	54	14920-0862	14720-0862
8	11.2	11.2	8.2	70	61	14920-0870	14720-0870
9	12.4	12.4	9.4	44	36	14920-0944	14720-0944
9	12.4	12.4	9.4	62	55	14920-0962	14720-0962
9	12.4	12.4	9.4	75	65	14920-0975	14720-0975
10	13.4	13.4	10.4	44	36	14920-1044	14720-1044
10	13.4	13.4	10.4	62	55	14920-1062	14720-1062
10	13.4	13.4	10.4	85	75	14920-1085	14720-1085
11	14.4	14.4	11.4	44	36	14920-1144	14720-1144
11	14.4	14.4	11.4	62	55	14920-1162	14720-1162
11	14.4	14.4	11.4	87	76	14920-1187	14720-1187
12	15.4	15.4	12.4	44	36	14920-1244	14720-1244
12	15.4	15.4	12.4	62	55	14920-1262	14720-1262
12	15.4	15.4	12.4	90	77	14920-1290	14720-1290
13	16.6	16.6	13.0	44	37	14920-1344	14720-1344
13	16.6	16.6	13.0	62	55	14920-1362	14720-1362
13	16.6	16.6	13.0	90	77	14920-1390	14720-1390

O.D. Tip = Outer diameter of the outer tube at the tip of the cannula · O.D. neck flange = Outer diameter behind the neck flange · I.D. = Inner diameter at the tip of the cannula  
 Length = Length over the outer curve · Length of mean radius = Length on midline, from flange to tip

STERILE EO

Sterilized using ethylene oxid



Keep away from sunlight



Single patient use



Consult instructions for use

Rx only

CAUTION:  
 United States Federal law restricts this device to sale, distribution and use by or on order of a physician or a licensed practitioner (in the US only).

CE marking identification number of the notified body

Product made in Germany