

HUMIDOFIX®

FILTER CASSETTE

- ▶ HME filter
- ▶ Filter the breathing air, thereby reducing intrusion of particles into the patient's airways
- ▶ Reduce the formation of viscous secretions in the lungs
- ▶ Store the moisture and warmth of the exhaled air in the filter medium and emit these to the inhaled air again during inspiration
- ▶ 15 mm standard connector fits to all tracheostomy tubes with 15 mm connector (UNI)
- ▶ Two paper filters integrated in one housing
- ▶ Low profile
- ▶ Adapter with oxygen supply tube available as separate accessory
- ▶ Individually packed in convenient peel-packs



Order information:

HUMIDOFIX®, sterile, pack of 30

REF 46440

Adapter for oxygen supply for HUMIDOFIX®, sterile, pack of 10

REF 46441