

HUMIDUAL®

sterile

REF 46470




- ▶ HME filter cassette
- ▶ 15 mm-Connector (UNI) - to connect to tracheostomy tubes with 15 mm connector
- ▶ Two open pore polyurethane foam filters, located in plastic housing in two opposite directions; double filtration capacity
- ▶ With oxygen supply port
- ▶ With overpressure valve to allow unimpaired coughing

SCOPE OF DELIVERY

- ▶ Instruction for use

HUMIDUAL®
1 x HME filter cassette

REF 46470

Product	Unit	REF
HUMIDUAL®	30 pieces, individually packed	46470

CE 0482