

LARYVOX® HANDS-FREE SET

Includes certain key devices for laryngectomized patients with voice prosthesis:

- ▶ LARYVOX® Hands-Free Valve Kombi HME Low Resistance
Adjustable speaking valve with 22 mm adapter (KOMBI) to connect an HME; fits into all LARYVOX® Tape base plates and LARYNGOTEC® silicone tubes and LARYNGOTEC® stoma buttons with 22 mm adapters (KOMBI)
- ▶ LARYVOX® Duo Brush Set with double-tip:
6 pieces (one box), size 10; for gentle cleaning of a voice prosthesis, with brush and sponge and LARYVOX® Clip to attach a stoma lamp; with mirror attached to the box
- ▶ LARYVOX® HME
5 pieces, low profile, inconspicuous, skin coloured, with stable ring housing for finger occlusion without pressure
- ▶ LARYVOX® Tape Comfort XL Oval
5 pieces, with very high adhesive strength for irregular tracheostomas and for physical activities
- ▶ OPTIGARD® Skin Protection Film
5 pieces, for preparative skin care when using base plates, helps to protect the skin against irritations and increases the adhesive power of the base plate
- ▶ OPTICLEAR® Adhesive Remover
5 pieces, to quickly and gently remove residual skin adhesives after the removal of the base plate

For all components of this kit separate and detailed data sheets are available

Please refer to the pictograms in the respective instructions of each product.

Order information:

LARYVOX® Hands-Free Set

REF 21986