


TRACHLINE® UNI LINGO-PHON

1 ICS + 1 ICU + 1 ICSU + SV

- ▶ Tracheostomy tube for long-term use, max. service life up to 29 days
- ▶ Made of thermosensitive medical synthetic materials
- ▶ Soft, pliable, lightweight
- ▶ Transparent, inconspicuous design
- ▶ Rounded tube tip to prevent irritation of the mucous membranes
- ▶ Bayonet locks enable simple, secure attachment of the inner cannula and accessories
- ▶ Set consists of a sieved inner cannula (ICS), a sieved inner cannula with 15 mm connector (ICSU), an inner cannula with 15 mm connector (ICU), a 15 mm connector (UNI), a cough cap, a sealing cap, a locking ring, a speaking valve, an adjustable tube holder, an obturator as well as a cannula passport
- ▶ TRACHLINE® PHON VALVE (speaking valve) and FAHL® UNI SEAL (sealing cap) are also available separately



Order information:

TRACHLINE® UNI LINGO-PHON
1 ICS + 1 ICU + 1 ICSU + SV

Size 5-10



REF 13804

FAHL® UNI SEAL

REF 13830

TRACHLINE® PHON VALVE

REF 27141

When ordering, please state the size after the item number!

Size	O.D. neck flange mm	I.D. O.C. mm	I.D. I.C. mm	Angel ca.	Length mm
5	9,1	6,5	5.0	115°	62
6	10,2	7,5	6.0	110°	70
7	11,2	8,5	7.0	105°	72
8	12,2	9,5	8.0	100°	80
9	13,1	10,5	9.0	95°	82
10	14,3	11,5	10.0	90°	84

O.D. neck flange = Outer diameter behind the neck flange · I.D. = Inner diameter at the tip of the cannula · O.C. = Outer cannula · I.C. = Inner cannula · Length = Length over middle curve