



# TRACHLINE® PHON VALVE

Speaking valve

- ▶ Very small, inconspicuous design and lightweight speaking valve
- ▶ Silicone valve
- ▶ The plastic housing is equipped on the patient side with a bayonet closure and thereby allows connection to a TRACHLINE® tracheostomy tube
- ▶ To be used by tracheotomized patients, not to be used by laryngectomized patients
- ▶ To be used only and exclusively in combination with a phonation cannula
- ▶ Individually packed



**Order information:**

TRACHLINE® PHON VALVE

Pack of 3

REF 27141