


# SENSOFOAM® / SENSOFOAM® PAD

- ▶ High-quality foam compress with high absorption capacity
- ▶ Made of a 3-layer foam material to prevent adhesion to the skin and to absorb secretions.
- ▶ Due to its reduced adhesion to the skin the compress may be removed easily
- ▶ Transmissible to air and humidity to prevent skin maceration and inflammation
- ▶ With slit to allow easy exchange of the compress while the tracheostomy tube remains in situ
- ▶ Stabilizes tight seal of the tracheostoma and a tight fit of the tracheostomy tube in the trachea
- ▶ Two variants: large and small
- ▶ SENSOFOAM® PAD with the shape of the flange of a tracheostomy tube to cushion the pressure of the neck-flange to reduce pain and the risk of pressure sores
- ▶ Sterile to allow immediate post-operative use



	Sterile
	Do not reuse
	Consult instructions for use
	CE marking

## Order information:

**SENSOFOAM® - sterile**, 10 x 9 x 0.5 cm, pack of 10

**REF 30860**

**SENSOFOAM® - sterile, small**, 6.5 x 6.5 x 0.5 cm, pack of 10

**REF 30861**

**SENSOFOAM® PAD - steril**, 7.8 x 3.8 cm, pack of 10

**REF 30862**