



- ▶ Five types of non-fenestrated inner cannulas for DURATWIX® tracheostomy tubes, each available in two lengths:
  1. With rotating 15 mm connector (ICV) to buffer tension and torque forces caused by the ventilation tubing (REF 19821 and REF 19822 short)
  2. With fixed 15 mm connector (ICU) for attachment of an HME (REF 19811 and REF 19812 short)
  3. With flat profile (IC) if no connection to any peripheral device (ventilation equipment, or HME) is desired (REF 19801 and REF 19802 short)
  4. With 22 mm adapter (ICK) to connect HMEs with 22 mm adapter (REF 19831 and REF 19832 short)
  5. Radiopaque, with fixed 15 mm connector (ICU) (REF 19803)
- ▶ Double colour coding:  
Connectors and/or connector base: white colour (turquoise colour for fenestrated tubes)  
Tube shaft: white colour for standard length variant; turquoise colour for short variant
- ▶ Size to be selected according to the size of the DURATWIX® outer cannula.

## Order information:

**DURATWIX® ICV**, SIZE 6-10, pack of 3

**DURATWIX® ICV short**, SIZE 6-10, pack of 3

**DURATWIX® ICU**, SIZE 6-10, pack of 3

**DURATWIX® ICU short**, SIZE 6-10, pack of 3

**DURATWIX® IC**, SIZE 7-10, pack of 3

**DURATWIX® IC short**, SIZE 7-10, pack of 3

**DURATWIX® ICK**, SIZE 7-10, pack of 3

**DURATWIX® ICK short**, SIZE 7-10, pack of 3

**DURATWIX® ICU contrast**, SIZE 6-10, pack of 3

	<b>REF 19821</b>
	<b>REF 19822</b>
	<b>REF 19811</b>
	<b>REF 19812</b>
	<b>REF 19801</b>
	<b>REF 19802</b>
	<b>REF 19831</b>
	<b>REF 19832</b>
	<b>REF 19803</b>

**STERILE EO** Sterilized using ethylene oxide

Do not re-sterilize

Keep away from sunlight

Single patient use

Consult instructions for use

**Rx only** CAUTION: United States Federal law restricts this device to sale, distribution and use by or on order of a physician or a licensed practitioner (in the US only).

**CE 0482** CE marking identification number of the notified body

**Product made in Germany**

When ordering, please state the size after the item number!