

COMBIPHON® SLIM

REF 27108



REF 27109



- ▶ Very small, inconspicuous design, transparent and lightweight speaking valve
- ▶ Silicone valve
- ▶ Fits to 15 mm-connector (UNI) devices
- ▶ To be used by tracheotomized patients, not to be used by laryngectomized patients
- ▶ To be used only and exclusively in combination with a phonation cannula
- ▶ **COMBIPHON® SLIM O₂** with oxygen supply port
- ▶ Individually packed

SCOPE OF DELIVERY

- ▶ Instruction for use
- | | |
|--------------------------------------|------------------|
| COMBIPHON® SLIM | REF 27108 |
| 1 x speaking valve | |
| or | |
| COMBIPHON® SLIM O₂ | REF 27109 |
| 1 x speaking valve | |

Product	Unit	REF
COMBIPHON® SLIM	1 piece	27108
COMBIPHON® SLIM O₂	1 piece	27109



UNI

O₂

The **COMBIPHON® SLIM** and **SLIM O₂** speaking valve must only be used by patients with partially or fully preserved larynx!