

NEBUJUNIOR® BEAR

INHALATION DEVICE FOR CHILDREN

- ▶ The NEBUJUNIOR® Bear is a modern inhalation device for compressed-air operated humidification of breathing air
- ▶ The fine aerosol mist enables optimum humidification of the respiratory tract
- ▶ The optical appearance, which appeals to children, reduces dislikes / fears associated with inhalation
- ▶ The NEBUJUNIOR® Bear can easily be disassembled into the individual components for maintenance and cleaning
- ▶ The NEBUJUNIOR® Bear inhalation device is supplied with a suitable storage bag in addition to the complete set of accessories (consisting of power supply unit, nebulizer, air tube, inhalation mask for adults, inhalation mask for children, air filter, mouthpiece and nose piece).
- ▶ Thus, inhalation is alternatively possible through the mouthpiece, the nose piece or the inhalation mask
- ▶ A special adapter tube with neck mask is available as special accessory for patients with tracheostoma
- ▶ Also available as penguin (NEBUJUNIOR® Penguin / REF 50400)



Technical Data

Power supply	AC 230 V, 50 Hz
Power consumption	80 VA +/- 20 VA
Max. filling quantity (medication)	13 ml
Particle size aerosol	0.5 - 10 µ m
Particle size acc. to MMAD	smaller than 4 µ m
Sound level	Ca. 60 dBA
Nebulization capacity	0.2 ml/min
Pressure range of compressor	29 - 43 Psi (2.0 - 3.0 bar)
Operating Pressure of compressor	8 - 16 Psi (0.7 - 1.1 bar)
Air flow	6 ~ 8 l/min
Operating temperature	10°C - 40°C (50°F - 104°F)
Rel. Air humidity during operation	10 - 95 % RH
Storage temperature	-20°C - 70°C (-4°F - 158°F)
Rel. Air humidity during storage	10 - 95 % RH

Technical Data

Dimensions	175 x 120 x 180 mm / (6.9" x 4.7" x 7.1")
Weight	1.6 kg (3.5 Lbs)
Device class	Type BF
IP class	II
Operation period	Short term operation – max. 30 minutes

Scope of delivery

1	Mains connection cable
1	Nebuliser
1	Air tube
1	Mouthpiece
1	Nosepiece
5	Filter
1	Inhalation mask for adults
1	Inhalation mask for children
1	Bag



Observe instructions for use!



Device of protection class II



Typ BF applied part



Please observe the laws and regulations which apply to the disposal of electrical and electronic devices in the respective countries.



CE marking identification number of the notified body



IP classification

Order information:

NEBUJUNIOR® BEAR
Tracheostoma mask adapter
LARYVOX® Mask

REF 50300

REF 50110

REF 75200