

# OPTISWAB® SET LARGE / OPTISWAB® SET XL

REF 31015



REF 31020



- ▶ Both sets for cleaning the tracheostomy tubes include:
  - **OPTIBRUSH® CONT XL** Cleaning tub (800 ml) (REF 31200) with large diameter to enable cleaning of the tracheostomy tubes. The cleaning tub contains a removable sieve insert which allows moving the tubes in the solution.
  - **OPTIBRUSH® CLEAN** cleaning powder (400 g) (REF 31100): for easy and accurate preparation of a cleaning solution; simply dissolve in lukewarm water. For thorough removal of adhering blood and secretion residues on the cannula.
- ▶ **OPTISWAB® SET LARGE (REF 31015):**
  - **OPTIBRUSH® SWAB** cleaning swabs (REF 31910) for the gentle cleaning of tracheostomy tubes. Flexible plastic swab with absorbent polyurethane foam cap. No fluff or fibres. This reduces the risk of tracheal irritation. Particularly suitable for tracheostomy tubes made of soft, sensitive plastic, especially for paediatric tracheostomy tubes.
- ▶ **OPTISWAB® SET XL (REF 31020):**
  - **OPTIBRUSH® SWAB XL** (REF 31920) cleaning swabs for the gentle cleaning of tracheostomy tubes. Flexible plastic swab with absorbent polyurethane foam cap. No fluff or fibres. This reduces the risk of tracheal irritation. Particularly suitable for tracheostomy tubes made of soft, sensitive plastic.

## SCOPE OF DELIVERY

- ▶ Instruction for use

### OPTISWAB® SET LARGE REF 31015

- 1 x OPTIBRUSH® CONT XL cleaning tub
- 1 x OPTIBRUSH® CLEAN cleaning powder (400 g)
- 4 x OPTIBRUSH® SWAB cleaning swab

or

### OPTISWAB® SET XL REF 31020

- 1 x OPTIBRUSH® CONT XL cleaning tub
- 1 x OPTIBRUSH® CLEAN cleaning powder (400 g)
- 4 x OPTIBRUSH® SWAB XL cleaning swab

All set components are also available individually.

## Product

OPTISWAB® SET LARGE

## REF

31015

OPTISWAB® SET XL

31020