




# LARYVOX<sup>®</sup> HME HIGHFLOW

- ▶ Special WIDE-PORE polyurethane foam to reduce respiration resistance
- ▶ Can hence be used up to 24 hours a day, even during physical activities
- ▶ Skin colour to provide inconspicuous design
- ▶ Lightweight; low profile
- ▶ Stable housing – can be closed by light finger occlusion to enable laryngectomized patients with a voice prosthesis to speak
- ▶ With 22 mm adapter (KOMBI)
- ▶ Pack size of 30 each



Keep away from sunlight



Keep dry



Do not reuse



Consult instructions for use



CE marking

## Order information:

LARYVOX<sup>®</sup> HME HighFlow

REF 49810